

# A STEP ABOVE THE REST

The only wireless, flat panel DR  
optimized for podiatry care



The Acuity DR 1012 Wireless DR is the next generation of digital radiology optimized specifically for podiatry care. While providing unmatched bone detail with industry-leading resolution, the Acuity DR 1012 provides a single solution to all podiatry examinations.



**Totally  
Wireless**



**No Cables**



**Long Battery  
Performance:  
630 images**



**Lightweight:  
4.1 lbs**



**Review HD  
processed  
image in as little  
as 6 seconds.**

**3x**

**Patient  
Throughput**



# No Cables, Any X-Ray, All Wireless

A reliable, yet simple Podiatry DR Solution

Acuity<sub>DR</sub>

Company A

# TRUVIEW ART

ADVANCED IMAGE RECONSTRUCTION TECHNOLOGY



AccuVue<sub>MED</sub>

10" x 12"



- 10 point touchscreen compatibility
- Powered by Radmedix's proprietary TruView Art image processing
- Fully customizable menus and exams
- DICOM Compliant
- Intelligent auto-advancing exams for patient care
- MWL, DICOM Send, Export, & CD Burn Included
- Built in patient x-ray and position guide

Leverage the Power of the AccuVue Cloud

AccuVue  
CLOUD



- *Computing and storage power over the internet.*
- *Unlimited storage and viewing without the cost to upgrade hardware.*
- *View and manipulate hundreds of large images using hand-held mobile devices like the iPad, iPhone and other tablets.*
- *Access images securely anytime, anywhere via cellular or Wi-Fi network access.*



## Included Guided Podiatry Tools

- *Hallux interphalangeus*
- *Hallux Abductus*
- *First Intermetatarsal*
- *P.A.S.A.*
- *D.A.S.A.*
- *Metatarsus Adductus*
- *Total Adductory Angle*
- *Metatarsal Protrusion*
- *Talo-Calcaneal*