

WBS-TM

Manual Tilting Wall Stand



Specification

Cassette stroke	1,400mm (400~1,800mm from floor to Bucky center)
Lowest tilting position	650mm from floor while tilting at 90 degrees
Bucky Type	Oscillating
Grid	FD 40~72inch (optimal FD 150cm), 103 lpi, ratio 10:1
Tilting Angle	-30 ~ +90 degrees
Inherent Filtration	0.5mmAl at 100kV
Lock(Brake)	EM Lock, Switch on/off (Vertical, Tilting)
Vertical Balance	Counter Weight
Tilting Balance	Spring
Electrical Rating	24VDC, 2.1A
Dimension / Weight	Max. 2,170(H) x 721.5(W) x 732(D) mm / 164kg(353lbs)
Option	Overhead handgrip