



# Sponges & Bolsters

2020  
PRODUCT  
GUIDE

## CHECK OUT WHAT'S NEW FOR 2020

Techno-Aide is ready to kick 2020 off with some amazing new products! Here are a few key highlights you can't afford to miss!

### **New Closed Cell Sponge Expansion | PG.3**

Techno-Aide now offers over 30 standard closed cell sponges in the most popular shapes and sizes! Shop our enormous in-stock inventory of individual closed cell sponges and kits!

#### **Closed Cell foam is the future of positioning sponges and the facts speak for themselves**


- Cleanable, non-porous surface
- Lighweight and comfortable
- Completely Artifact Free
- 100% radiolucent foam
- Extremely durable




### **Completely Sealed Vinyl Bolsters | PG.17**

Our vinyl fabric is now sewn and welded at all seams to offer a completely sealed and secure bolster! These premium sealed corners inhibit the growth of any harmful bacteria inside and out. Compared to other unsealed and unsanitary, "hook and loop" closures, these bolsters are the clear and easy choice.

Completely Sealed



Unsealed & Unsanitary



## Picking the Perfect Positioning Aid

Techno-Aide is here to help you find the perfect positioning aid! Our sponges are specially designed to meet any need! Below will be everything you need to know about our sponge options and how to choose the perfect fit.

### STEP 1: WHAT TYPE OF SPONGE DO YOU NEED?


#### Uncoated Foam

Basic (core), open-cell foam. Non-coated, non-cleanable surface. The least durable and least expensive foam.



#### Coated Foam

Fluid-resistant, cleanable coating over standard (core) foam core. Minimal shadowing with low attenuation. Light green.




#### Closed Cell Foam

**Our most popular foam!** Artifact free, cleanable surface, lightweight and durable. Our top option! Charcoal grey.



#### Vinyl Covered Foam

Firm foam with a heavy-duty vinyl covering. Ideal for extreme patient comfort and durability. Blue vinyl.



Go to steps 2 & 3!

Skip to step 4!

### STEP 2: FOAM TYPES

#### Standard Foam

Standard open-cell foam for general positioning. Light and comfortable option for basic imaging.



#### Firm/Bariatric Foam


A higher density open-cell foam than standard foam. Typically used to support larger/heavier/bariatric patients for positioning.



### STEP 3: ARTIFACTS?

#### Non-Flared Sides

**Non-flared sponges may show artifacts under imaging, specifically when coating is applied.**



#### Flared Sides

**Stealth™** sponges are designed with proprietary angles to offer **artifact-free** viewing. Flared sides will increase size!



### STEP 4: PICK A SHAPE AND SIZE!

Different imaging procedures require different designs. Techno-Aide is proud to offer the largest selection of standard and specialty sponges available. We offer hundreds of sponges designed for spinal, epidural, cardiac, echocardiography, extremity, cervical, AP/PA, decubitus, ultrasound, phlebotomy, CT, podiatry and many more! Starting on page 3, shop the rest of this catalog for our most popular options or submit for a custom quote today!

## NEW Closed Cell Positioning Sponges

Closed cell foam is the future of radiology. This cleanable, non-porous foam is completely radiolucent. The charcoal grey, low density foam is extremely durable and comfortable for patients. Closed cell sponges offer 100% artifact free imaging.

### STOCKED SPONGES

| CAT.# | Description       | Size (W x L x H)     | Price    |
|-------|-------------------|----------------------|----------|
| YXBC  | 45° Wedge         | 7" x 7" x 7"         | \$89.00  |
| YXBB  | 45° Wedge         | 7" x 14" x 7"        | \$102.00 |
| YXBD  | 45° Wedge         | 7" x 21" x 7"        | \$85.00  |
| YXEA  | 45° Wedge         | 9.5" x 24" x 9.5"    | \$153.00 |
| YX16  | 45° Wedge         | 9" x 30" x 9"        | \$137.00 |
| YXBU  | 30°/60°/90° Wedge | 24" x 10" x 6"       | \$132.00 |
| YXBE  | 34° Wedge         | 7" x 4" x 3"         | \$139.00 |
| YXBJ  | 12° Wedge Set     | 9.5" x 7.25" x 2.25" | \$167.00 |
| YXBI  | 15° Wedge Set     | 7" x 10" x 2"        | \$114.00 |
| YXAJ  | Knee Wedge        | 22" x 13.5" x 7"     | \$191.00 |
| YXAK  | 7" Head Circle    | 7" Diameter          | \$15.00  |
| YXAA  | Oblique Finger    | 5" x 5" x 3"         | \$20.00  |
| YXAT  | Angular Bolster   | 16" x 12" x 7"       | \$146.00 |
| YXAY  | Angular Bolster   | 16" x 12.5" x 7.25"  | \$169.00 |
| YXUZ  | Ultrasound Wedge  | 20" x 20" x 6"       | \$175.00 |
| YXBV  | Rectangle         | 2" x 6" x 6"         | \$20.00  |
| YXMB  | Rectangle         | 1" x 8" x 10"        | \$24.00  |
| YXBN  | Rectangle         | 3" x 14" x 18"       | \$47.00  |
| YXBL  | Rectangle         | 3" x 7.5" x 10"      | \$79.00  |
| YXGM  | Rectangle         | 4" x 18" x 24"       | \$72.00  |

**Stocked sponges  
are marked Blue!**

### NON-STOCKED SPONGES

| CAT.# | Description   | Size (W x L x H)         | Price    |
|-------|---------------|--------------------------|----------|
| YX278 | Square        | 12" x 12" x 12"          | \$109.00 |
| YX277 | Square        | 10" x 10" x 10"          | \$99.00  |
| YX100 | Rectangle     | 10" x 14" x 17"          | \$128.00 |
| YX169 | Rectangle     | 6" x 18" x 42"           | \$249.00 |
| YX73  | Rectangle     | 6" x 18" x 22"           | \$153.00 |
| YX126 | Rectangle     | 6" x 8" x 10"            | \$38.00  |
| YX105 | Rectangle     | 6" x 8" x 12"            | \$40.00  |
| YX70  | Rectangle     | 5" x 10" x 24"           | \$66.00  |
| YX58  | Rectangle     | 5" x 14" x 17"           | \$75.00  |
| YX67  | Rectangle     | 5" x 3" x 7"             | \$42.00  |
| YX71  | Rectangle     | 4" x 10" x 20"           | \$77.00  |
| YX62  | Rectangle     | 4" x 10" x 14"           | \$65.00  |
| YX03  | Rectangle     | 4" x 6" x 6"             | \$29.00  |
| YX250 | Rectangle     | 3" x 9" x 24"            | \$50.00  |
| YX204 | Rectangle     | 3" x 12" x 14"           | \$40.00  |
| YX156 | Rectangle     | 3" x 6" x 6"             | \$24.00  |
| YX266 | Rectangle     | 4" x 18" x 36"           | \$140.00 |
| YX06  | Rectangle     | 2" x 8" x 24"            | \$34.00  |
| YX246 | Rectangle     | 2" x 14" x 18"           | \$51.00  |
| YX29  | Rectangle     | 2" x 8" x 16"            | \$29.00  |
| YXGF  | Rectangle     | 2" x 10" x 12"           | \$26.00  |
| YXBT  | Rectangle     | 2" x 4" x 6"             | \$19.00  |
| YXBR  | Rectangle     | 1" x 4" x 6"             | \$18.00  |
| YXER  | Rectangle     | 1" x 4" x 12"            | \$19.00  |
| YXBQ  | Rectangle     | 0.5" x 4" x 6"           | \$19.00  |
| YXBH  | 16° Wedge Set | 3.5" x 8" x 7"           | \$164.00 |
| YXBG  | 21° Wedge     | 4.625" x 10" x 10"       | \$94.00  |
| YXBF  | 27° Wedge     | 6" x 10" x 10"           | \$99.00  |
| YX51  | 33° Wedge     | 6.6" x 9" x 12"          | \$112.00 |
| YX167 | 45° Wedge     | 10" x 10" x 14"          | \$135.00 |
| YX56  | 45° Wedge     | 9.5" x 13.75" x 24"      | \$180.00 |
| YX209 | 45° Wedge     | 6" x 12" x 38"           | \$256.00 |
| YX75  | 45° Wedge     | 10" x 10" x 16"          | \$152.00 |
| YX76  | 45° Wedge     | 3" x 8" x 3"             | \$65.00  |
| YXCA  | Trough Wedge  | Small - 4" x 10" x 3"    | \$75.00  |
| YXCB  | Trough Wedge  | Medium - 7" x 15" x 5"   | \$90.00  |
| YXCC  | Trough Wedge  | Large - 11" x 24" x 7"   | \$150.00 |
| YXCD  | Trough Wedge  | X-Large - 14" x 36" x 9" | \$276.00 |

### NON-STOCKED SPONGES

| CAT.# | Description               | Size (W x L x H)   | Price    |
|-------|---------------------------|--------------------|----------|
| YXAB  | Finger                    | 11" x 5" x 3"      | \$113.00 |
| YXAC  | Finger                    | 7.5" x 4" x 1.875" | \$87.00  |
| YXAD  | Finger                    | 13" x 5.5" x 2.75" | \$75.00  |
| YXAE  | Finger                    | 7.25" x 3" x 2.75" | \$68.00  |
| YX126 | Knee Block                | 6" x 8" x 10"      | \$38.00  |
| YXAG  | Bilateral Leg Immobilizer | 6" x 14" x 16"     | \$149.00 |

## NEW Orthopedic Leg Length Discrepancy Block Kit

Convenient and practical method for the assessing and testing of patients with leg length discrepancy. Nine closed cell foam blocks are graded by half centimeters (1cm to 5cm). The block is placed under the foot of the shorter leg, to achieve proper position of the pelvis for leg, spine, and pelvis x-rays in the standing position. The image is best performed with blocks under the short limb to level the pelvis, to avoid unwanted adjustments by the patient to their hips, knees or ankles.

YZLL | \$445.00

| CAT#   | Width x Length (in) | Height (cm) | QTY |
|--------|---------------------|-------------|-----|
| YXLL01 | 10" x 14"           | 1 cm        | 1   |
| YXLL15 | 10" x 14"           | 1.5 cm      | 1   |
| YXLL02 | 10" x 14"           | 2 cm        | 1   |
| YXLL25 | 10" x 14"           | 2.5 cm      | 1   |
| YXLL03 | 10" x 14"           | 3 cm        | 1   |
| YXLL35 | 10" x 14"           | 3.5 cm      | 1   |
| YXLL04 | 10" x 14"           | 4 cm        | 1   |
| YXLL45 | 10" x 14"           | 4.5 cm      | 1   |
| YXLL05 | 10" x 14"           | 5 cm        | 1   |




## Podiatry Axial/Sesamoid Closed Cell Sponge

Designed specifically to view metatarsals, sesamoids and plantar condyles. Made of our 2 lb., non-porous, cleanable, fluid resistant closed cell foam. Charcoal grey in color.

#YXNZ | \$192.00

8" x 13" x 4"



## NEW Closed Cell Specialty Knee Bolster

Designed to support patients upper thigh and knees during imaging procedures. Helps to alleviate pressure on the lower back providing extra comfort.

#YXKB | \$233.00

16" x 25" x 8.25"



## NEW Weight Bearing Receptor Holder

**NEW** Two options to accommodate nearly EVERY receptor type!  
Economic solution for weight bearing examinations. Durable 6 lb closed cell provides support for weight bearing studies.

| CAT#          | Sponge Dimensions                | Typical Receptor Size | Price           |
|---------------|----------------------------------|-----------------------|-----------------|
| YXWB          | 11.5"(W) x 17.5"(L) x 4"(H)      | 10" x 12"             | \$146.00        |
| <b>YXW-18</b> | <b>18"(W) x 17.5"(L) x 4"(H)</b> | <b>14" x 17"</b>      | <b>\$210.00</b> |



**Channel is 1 3/16" wide**  
to accommodate most  
imaging receptors.

## NEW Closed Cell Positioning Sponge Kits

Closed cell foam is the future of radiology. This cleanable, non-porous foam is completely radiolucent. The charcoal grey, low density foam is extremely durable and comfortable for patients. Closed cell sponges offer 100% artifact free imaging.

### STOCKED KITS

#### GENERAL KIT 1 | YZCC - \$199.00

| CAT.# | Description           | Size (W x L x H) | QTY |
|-------|-----------------------|------------------|-----|
| YXBD  | 45° Degree Body Wedge | 7" x 21" x 7"    | 1   |
| YXAA  | Oblique Finger Block  | 5" x 5" x 3"     | 1   |
| YXAK  | Circular Head Block   | 7" Diam.         | 1   |
| YXBI  | 15° Degree Wedge Set  | 7" x 10" x 2"    | 1   |
| YXBV  | Rectangle Block       | 6" x 6" x 2"     | 1   |
| YXMB  | Rectangle Block       | 8" x 10" x 2"    | 1   |
| YXBN  | Rectangle Block       | 14" x 18" x 3"   | 1   |

#### GENERAL KIT A | YZGA - \$340.00

| CAT.# | Description           | Size (W x L x H) | QTY |
|-------|-----------------------|------------------|-----|
| YXBC  | 45° Degree Body Wedge | 7" x 7" x 7"     | 1   |
| YXAK  | Circular Head Block   | 7" Diam.         | 1   |
| YXBI  | 15° Degree Wedge Set  | 7" x 10" x 2"    | 1   |
| YXBE  | 34° Degree Wedge Set  | 7" x 4" x 3"     | 1   |
| YXBL  | Rectangle             | 4" x 18" x 24"   | 1   |

#### GENERAL KIT B | YZGB - \$536.00

| CAT.# | Description           | Size (W x L x H)  | QTY |
|-------|-----------------------|-------------------|-----|
| YXBC  | 45° Degree Body Wedge | 7" x 7" x 7"      | 1   |
| YXAA  | Oblique Finger Block  | 5" x 5" x 3"      | 1   |
| YXAK  | Circular Head Block   | 7" Diam.          | 1   |
| YXBI  | 15° Degree Wedge Set  | 7" x 10" x 2"     | 1   |
| YXBE  | 34° Degree Wedge Set  | 7" x 4" x 3"      | 1   |
| YXBJ  | 12° Degree Wedge Set  | 9.5" x 7.25" x 3" | 1   |
| YXBL  | Rectangle Block       | 3" x 7.5" x 10"   | 1   |

### NON-STOCKED KITS

#### GENERAL KIT 2 | YZG2 - \$715.00

| CAT.# | Description           | Size (W x L x H) | QTY |
|-------|-----------------------|------------------|-----|
| YXAK  | Circular Head Block   | 7" Diam.         | 1   |
| YXMB  | Rectangle Block       | 8" x 10" x 2"    | 2   |
| YXBE  | 34° Degree Wedge Set  | 7" x 4" x 3"     | 1   |
| YXBH  | 16° Degree Wedge Set  | 3.5" x 8" x 7"   | 1   |
| YXBC  | 45° Degree Body Wedge | 7" x 7" x 7"     | 1   |

#### TRIANGLE WEDGE KIT | YZTK - \$685.00

| CAT.# | Description        | Size (W x L x H) | QTY |
|-------|--------------------|------------------|-----|
| YXBG  | 21° Triangle Wedge | 4.6" x 10" x 10" | 2   |
| YX51  | 33° Triangle Wedge | 6.6" x 9" x 12"  | 2   |
| YX167 | 45° Triangle Wedge | 10" x 10" x 14"  | 1   |
| YXBF  | 27° Triangle Wedge | 6" x 10" x 10"   | 2   |

**Stocked sponges  
are marked Blue!**

## NEW Closed Cell Positioning Sponge Kits

Closed cell foam is the future of radiology. This cleanable, non-porous foam is completely radiolucent. The charcoal grey, low density foam is extremely durable and comfortable for patients. Closed cell sponges offer 100% artifact free imaging.

### STOCKED KITS

#### CLINIC KIT A | YZCA - \$181.00

| CAT.# | Description           | Size (W x L x H) | QTY |
|-------|-----------------------|------------------|-----|
| YXBB  | 45° Degree Body Wedge | 7" x 14" x 7"    | 1   |
| YXAA  | Oblique Finger Block  | 5" x 5" x 3"     | 1   |
| YXAK  | Circular Head Block   | 7" Diam.         | 1   |
| YXBI  | 15° Degree Wedge Set  | 7" x 10" x 2"    | 1   |

#### CLINIC KIT B | YZCB - \$679.00

| CAT.# | Description          | Size (W x L x H)  | QTY |
|-------|----------------------|-------------------|-----|
| YXBD  | 45° Body Wedge       | 7" x 21" x 7"     | 1   |
| YXAA  | Oblique Finger Block | 5" x 5" x 3"      | 1   |
| YXAK  | Circular Head Block  | 7" Diam.          | 1   |
| YXBI  | 15° Wedge Set        | 7" x 10" x 2"     | 1   |
| YXBL  | Rectangle            | 3" x 7.5" x 10"   | 1   |
| YXBC  | 45° Body Wedge       | 7" x 7" x 7"      | 1   |
| YXBE  | 34° Wedge Set        | 7" x 4" x 3"      | 1   |
| YXBJ  | 12° Wedge Set        | 9.5" x 7.25" x 3" | 1   |

#### GENERAL KIT P | YZGP - \$148.00

| CAT.# | Description          | Size (W x L x H) | QTY |
|-------|----------------------|------------------|-----|
| YXBD  | 45° Body Wedge       | 7" x 21" x 7"    | 1   |
| YXAA  | Oblique Finger Block | 5" x 5" x 3"     | 1   |
| YXAK  | Circular Head Block  | 7" Diam.         | 1   |
| YXBI  | 15° Wedge Set        | 7" x 10" x 2"    | 1   |



#YXAK



#YXBN



#YXBI

### NON-STOCKED KITS

#### ORTHO KIT | YZOR - \$820.00

| CAT.# | Description          | Size (W x L x H)     | QTY |
|-------|----------------------|----------------------|-----|
| YXNZ  | Podiatry Axial Spoge | 8" x 13" x 4"        | 1   |
| YXEA  | 45° Wedge            | 9.5" x 24" x 9.5"    | 2   |
| YXAA  | Oblique Finger       | 5" x 5" x 3"         | 1   |
| YXBG  | 21° Wedge            | 4.6" x 10" x 10"     | 1   |
| YXAK  | 7" Head Circle       | 7" Diam.             | 1   |
| YXBJ  | Knee Wedge Set       | 9.5" x 7.25" x 2.25" | 1   |
| YXFB  | Semi-Round Cylinder  | 8" x 16" x 6"        | 1   |
| YX59  | Cylinder             | 24" x 8"             | 1   |

#### GENERAL KIT | YZVE - \$200.00

| CAT.# | Description        | Size (W x L x H) | QTY |
|-------|--------------------|------------------|-----|
| YXGM  | Rectangle          | 4" x 18" x 24"   | 1   |
| YX53  | 33° Triangle Wedge | 9.5" x 7.5" x 2" | 2   |
| YXBE  | 34° Wedge Set      | 7" x 4" x 3"     | 2   |
| YXBH  | 16° Wedge Set      | 3.5" x 8" x 7"   | 2   |

#### TROUGH KIT | YZ XV - \$592.00

| CAT.# | Description  | Size (W x L x H)         | QTY |
|-------|--------------|--------------------------|-----|
| YXCA  | Trough Wedge | Small - 4" x 10" x 3"    | 1   |
| YXCB  | Trough Wedge | Medium - 7" x 15" x 5"   | 1   |
| YXCC  | Trough Wedge | Large - 11" x 24" x 7"   | 1   |
| YXCD  | Trough Wedge | X-Large - 14" x 36" x 9" | 1   |

## Stealth Flared Positioning Sponges

Revolutionary design minimizes ghosting and line artifacts. Stealth positioning sponges are available coated and non-coated. May be larger than traditional style. Dimensions listed as "Width x Length x Height."

**Stocked sponges are marked Blue!**

### COATED SPONGES

#### STEALTH RECTANGLES

| CAT.# | Base (W x L x H)    | Top (W x L x H) | Price    |
|-------|---------------------|-----------------|----------|
| YCBL  | 10.25" x 12.5" x 3" | 7.5" x 10"      | \$63.00  |
| YCBN  | 16.25" x 20.5" x 3" | 14" x 18"       | \$161.00 |
| YCBO  | 13.5" x 13.5" x 4"  | 10" x 10"       | \$83.00  |
| YCBY  | 21.5" x 27.5" x 4"  | 18" x 24"       | \$233.00 |

#### STEALTH WEDGE SETS (TWO EACH)

| CAT.# | Degree | Top (W x L x H)       | QTY | Price    |
|-------|--------|-----------------------|-----|----------|
| YCBE  | 34°    | 5.5" x 9.75" x 3"     | 2   | \$55.00  |
| YCBJ  | 12°    | 10.5" x 9.25" x 2.25" | 2   | \$92.00  |
| YCBH  | 16°    | 8.5" x 11" x 3.5"     | 2   | \$108.00 |
| YCBI  | 15°    | 12.5" x 10" x 3"      | 2   | \$124.00 |

#### STEALTH WEDGES

| CAT.# | Degree      | Top (W x L x H)        | Price    |
|-------|-------------|------------------------|----------|
| YCBC  | 45°         | 10" x 13.5" x 7"       | \$65.00  |
| YCBF  | 27°         | 12.75" x 15.5" x 6"    | \$76.00  |
| YCBG  | 21°         | 12.25" x 14.25" x 4.5" | \$88.00  |
| YCBB  | 45°         | 10" x 20.25" x 7"      | \$101.00 |
| YCEB  | 45°         | 10" x 30" x 7.5"       | \$138.00 |
| YCEA  | 45°         | 13.75" x 33" x 9.5"    | \$135.00 |
| YCBD  | 45°         | 10" x 27.75" x 7"      | \$142.00 |
| YCEE  | 45°         | 13" x 38.25" x 9"      | \$145.00 |
| YCBU  | 30°/60°/90° | 10" x 29.5" x 6"       | \$196.00 |

### NON-COATED SPONGES

#### STEALTH RECTANGLES


| CAT.# | Base (W x L x H)    | Top (W x L x H) | Price   |
|-------|---------------------|-----------------|---------|
| YFBL  | 10.25" x 12.5" x 3" | 7.5" x 10"      | \$14.00 |
| YFBN  | 16.25" x 20.5" x 3" | 14" x 18"       | \$37.00 |
| YFBO  | 13.5" x 13.5" x 4"  | 10" x 10"       | \$26.00 |
| YFBY  | 21.5" x 27.5" x 4"  | 18" x 24"       | \$64.00 |

#### FLARED WEDGE SETS

| CAT.# | Degree | Top (W x L x H)       | QTY | Price   |
|-------|--------|-----------------------|-----|---------|
| YFBE  | 34°    | 5.5" x 9.75" x 3"     | 2   | \$15.00 |
| YFBJ  | 12°    | 10.5" x 9.25" x 2.25" | 2   | \$16.00 |
| YFBH  | 16°    | 8.5" x 11" x 3.5"     | 2   | \$21.00 |
| YFBI  | 15°    | 12.5" x 10" x 3"      | 2   | \$20.00 |

#### STEALTH WEDGES

| CAT.# | Degree      | Top (W x L x H)        | Price   |
|-------|-------------|------------------------|---------|
| YFBC  | 45°         | 10" x 13.5" x 7"       | \$21.00 |
| YFBF  | 27°         | 12.75" x 15.5" x 6"    | \$23.00 |
| YFBG  | 21°         | 12.25" x 14.25" x 4.5" | \$18.00 |
| YFBB  | 45°         | 10" x 20.25" x 7"      | \$26.00 |
| YFEB  | 45°         | 10" x 30" x 7.5"       | \$50.00 |
| YFEA  | 45°         | 13.75" x 33" x 9.5"    | \$44.00 |
| YFBD  | 45°         | 10" x 27.75" x 7"      | \$42.00 |
| YFEE  | 45°         | 13" x 38.25" x 9"      | \$43.00 |
| YFBU  | 30°/60°/90° | 10" x 29.5" x 6"       | \$33.00 |





## Stealth Flared Positioning Sponges

Revolutionary design minimizes ghosting and line artifacts. Stealth positioning sponges are available coated and non-coated. May be larger than traditional style. Dimensions listed as "Width x Length x Height."

**Stocked sponges  
are marked Blue!**

### COATED SPONGES

STEALTH FINGER BLOCKS

| CAT.# | Description/Use                          | Price   |
|-------|--|---------|
| YCAC  | Bilateral Pediatric Oblique Finger Block | \$32.00 |
| YCAA  | Adult Oblique Finger Block               | \$34.00 |
| YCAE  | Pediatric Lateral Finger Block           | \$45.00 |
| Y CAB | Bilateral Adult Oblique Finger Block     | \$59.00 |
| YCAD  | Adult Lateral Finger Block               | \$78.00 |

STEALTH EXTREMITY SPONGES

| CAT.# | Description            | Size (W x L x H)    | Price    |
|-------|------------------------|---------------------|----------|
| YCCB  | Medium Extremity Block | 11.5" x 19.75" x 5" | \$123.00 |

STEALTH KNEE SPONGES

| CAT.# | Description           | Size (W x L x H) | Price    |
|-------|-----------------------|------------------|----------|
| YCAH  | Knee Arthrogram Block | 14.5" x 19" x 8" | \$135.00 |
| YCAJ  | Knee Wedge Support    | 22" x 20" x 7"   | \$315.00 |

### NON-COATED SPONGES

STEALTH FINGER BLOCKS


| CAT.# | Description/Use                          | Price   |
|-------|--|---------|
| YFAC  | Bilateral Pediatric Oblique Finger Block | \$14.00 |
| YFAA  | Adult Oblique Finger Block               | \$14.00 |
| YFAE  | Pediatric Lateral Finger Block           | \$25.00 |
| YFAB  | Bilateral Adult Oblique Finger Block     | \$27.00 |
| YFAD  | Adult Lateral Finger Block               | \$28.00 |

STEALTH EXTREMITY SPONGES

| CAT.# | Description            | Size (W x L x H)    | Price   |
|-------|------------------------|---------------------|---------|
| YFCB  | Medium Extremity Block | 11.5" x 19.75" x 5" | \$42.00 |

STEALTH KNEE SPONGES


| CAT.# | Description           | Size (W x L x H) | Price    |
|-------|-----------------------|------------------|----------|
| YFAH  | Knee Arthrogram Block | 14.5" x 19" x 8" | \$48.00  |
| YFAJ  | Knee Wedge Support    | 22" x 20" x 7"   | \$109.00 |




#YCAD



#YCAH



#YCCB



#YCAJ

## Stealth Flared Positioning Sponges

Revolutionary design minimizes ghosting and line artifacts. Stealth positioning sponges are available coated and non-coated. May be larger than traditional style. Dimensions listed as "Width x Length x Height."

**Stocked sponges are marked Blue!**

### COATED SPONGES

#### STEALTH HEAD & NECK SPONGES

| CAT.# | Description             | Size (W x L x H)      | Price    |
|-------|-------------------------|-----------------------|----------|
| YCAK  | Circular Head Block     | 7-9" Diam. x 1.5" (H) | \$37.00  |
| YCAR  | Cervical Head Rest      | 7.5" x 16" x 4.5"     | \$83.00  |
| YCAL  | Head & Neck Support     | 14.5" x 16" x 6.5"    | \$88.00  |
| YCAN  | AP/PA Child Head Holder | 9.75" x 11" x 5"      | \$91.00  |
| YCAQ  | Child Head Immobilizer  | 6.75" x 11.5" x 5.25" | \$99.00  |
| YCAM  | AP/PA Adult Head Holder | 14" x 15" x 6.75"     | \$145.00 |
| YCAP  | Adult Head Immobilizer  | 9" x 15.25" x 7"      | \$176.00 |

#### STEALTH BARIATRIC RECTANGLE/DECUBITUS

| CAT.# | Base (W x L x H)   | Top (W x L x H) | Price    |
|-------|--------------------|-----------------|----------|
| YCBM  | 10.5" x 14.5" x 3" | 8" x 12"        | \$68.00  |
| YCBP  | 21.5" x 27.5" x 4" | 18" x 24"       | \$217.00 |

#### STEALTH BARIATRIC TORSO SPONGES

| CAT.# | Description             | Size (W x L x H)   | Price    |
|-------|-------------------------|--------------------|----------|
| YCAT  | Bariatric Torso Bolster | 18.5" x 16.5" x 7" | \$165.00 |

#### STEALTH BARIATRIC WEDGES

| CAT.# | Degree      | Top (W x L x H)  | Price    |
|-------|-------------|------------------|----------|
| YCBA  | 45°         | 10" x 27.5" x 7" | \$150.00 |
| YCEU  | 30°/60°/90° | 10" x 29.5" x 6" | \$218.00 |

### NON-COATED SPONGES

#### STEALTH HEAD & NECK SPONGES

| CAT.# | Description             | Size (W x L x H)      | Price   |
|-------|-------------------------|-----------------------|---------|
| YFAK  | Circular Head Block     | 7-9" Diam. x 1.5" (H) | \$15.00 |
| YFAR  | Cervical Head Rest      | 7.5" x 16" x 4.5"     | \$23.00 |
| YFAL  | Head & Neck Support     | 14.5" x 16" x 6.5"    | \$33.00 |
| YFAN  | AP/PA Child Head Holder | 9.75" x 11" x 5"      | \$29.00 |
| YFAQ  | Child Head Immobilizer  | 6.75" x 11.5" x 5.25" | \$37.00 |
| YFAM  | AP/PA Adult Head Holder | 14" x 15" x 6.75"     | \$47.00 |
| YFAP  | Adult Head Immobilizer  | 9" x 15.25" x 7"      | \$56.00 |

#### STEALTH BARIATRIC RECTANGLE/DECUBITUS

| CAT.# | Base (W x L x H)   | Top (W x L x H) | Price   |
|-------|--------------------|-----------------|---------|
| YBBM  | 10.5" x 14.5" x 3" | 8" x 12"        | \$24.00 |
| YBBP  | 21.5" x 27.5" x 4" | 18" x 24"       | \$96.00 |

#### STEALTH BARIATRIC TORSO SPONGES

| CAT.# | Description             | Size (W x L x H)   | Price   |
|-------|-------------------------|--------------------|---------|
| YBAT  | Bariatric Torso Bolster | 18.5" x 16.5" x 7" | \$59.00 |

#### STEALTH BARIATRIC WEDGES

| CAT.# | Degree      | Top (W x L x H)  | Price   |
|-------|-------------|------------------|---------|
| YBBA  | 45°         | 10" x 27.5" x 7" | \$46.00 |
| YFEU  | 30°/60°/90° | 10" x 29.5" x 6" | \$47.00 |



#YCAK



#YCAR



#YCAN



#YCAL



#YCAQ

## Non-Stealth Positioning Sponges

Non-Stealth positioning sponges are available coated and non-coated. Non-flared sponges may show line artifacts when coating is applied. Dimensions listed as "Width x Length x Height."

**Stocked sponges are marked Blue!**

### COATED SPONGES

### NON-COATED SPONGES

#### NON-STEALTH RECTANGLES

| CAT.# | Overall Size (W x L x H) | Price    |
|-------|--------------------------|----------|
| YCBR  | 4" x 6" x 1"             | \$14.00  |
| YCBS  | 4" x 6" x 1.5"           | \$24.00  |
| YCBT  | 4" x 6" x 2"             | \$30.00  |
| YCBV  | 6" x 6" x 2"             | \$36.00  |
| YCGG  | 8" x 18" x 1.5"          | \$91.00  |
| YCEO  | 14" x 21" x 1.5"         | \$130.00 |
| YCEM  | 15" x 24" x 2.5"         | \$201.00 |
| YCEH  | 5.5" x 20" x 4"          | \$68.00  |

#### NON-STEALTH RECTANGLES

| CAT.# | Overall Size (W x L x H) | Price   |
|-------|--------------------------|---------|
| YFBR  | 4" x 6" x 1"             | \$5.00  |
| YFBS  | 4" x 6" x 1.5"           | \$6.00  |
| YFBT  | 4" x 6" x 2"             | \$7.00  |
| YFBV  | 6" x 6" x 2"             | \$19.00 |
| YFGG  | 8" x 18" x 1.5"          | \$15.00 |
| YFEO  | 14" x 21" x 1.5"         | \$45.00 |
| YFEM  | 15" x 24" x 2.5"         | \$31.00 |
| YFEH  | 5.5" x 20" x 4"          | \$25.00 |

#### NON-STEALTH WEDGES

| CAT.# | Degree | Type     | Size (W x L x H)  | Price    |
|-------|--------|----------|-------------------|----------|
| YCFD  | 45°    | Spinal   | 7" x 21" x 7"     | \$120.00 |
| YCEG  | 65°    |          | 16" x 25" x 8.25" | \$195.00 |
| YCCL  | 21°    | Epidural | 24" x 27" x 7"    | \$217.00 |
| YCDB  | 30°    | Cardiac  | 28" x 27" x 16"   | \$264.00 |
| YCDA  | 45°    | Cardiac  | 19" x 28" x 14"   | \$289.00 |

#### NON-STEALTH WEDGES


| CAT.# | Degree | Type     | Size (W x L x H)  | Price    |
|-------|--------|----------|-------------------|----------|
| YFFD  | 45°    | Spinal   | 7" x 21" x 7"     | \$37.00  |
| YFEG  | 65°    |          | 16" x 25" x 8.25" | \$68.00  |
| YFCL  | 21°    | Epidural | 24" x 27" x 7"    | \$88.00  |
| YFDB  | 30°    | Cardiac  | 28" x 27" x 16"   | \$152.00 |
| YFDA  | 45°    | Cardiac  | 19" x 28" x 14"   | \$192.00 |

#### NON-STEALTH BARIATRIC WEDGES

| CAT.# | Degree | Type   | Size (W x L x H) | Price    |
|-------|--------|--------|------------------|----------|
| YCCG  | 45°    |        | 9" x 30" x 9"    | \$153.00 |
| YCHB  |        | Spinal | 7" x 24" x 7"    | \$182.00 |

#### NON-STEALTH BARIATRIC WEDGES


| CAT.# | Degree | Type   | Size (W x L x H) | Price   |
|-------|--------|--------|------------------|---------|
| YFCG  | 45°    |        | 9" x 30" x 9"    | \$47.00 |
| YFHB  |        | Spinal | 7" x 24" x 7"    | \$55.00 |



## Non-Stealth Specialty Traditional Sponges

Specialty shaped sponges may show line artifacts when coating is applied.

**Stocked sponges are marked Blue!**




COATED

NON-COATED

#YCUZ  
\$275.00

#YBUZ  
\$112.00

**ENDO-ULTRASOUND WEDGE**  
20 x 20 x 6



#YCEQ  
\$178.00

#YFEQ  
\$78.00

**PERINEAL WEDGE**  
20 x 20 x 3



#YCAV  
\$172.00

#YFAY  
\$56.00

**ANGULAR BOLSTER**  
14 x 22.5 x 7.25




## Non-Stealth Specialty Sponges

Specialty shaped sponges may show line artifacts when coating is applied.




#YCEZ  
\$102.00  
**PHLEBOTOMY WEDGE**  
5.5 x 25 x 5.5



#YCBW  
\$249.00  
**CT HEADREST**  
12.5 x 12 x 9.75  
(Neck Cradle 3.75" h)



#YCMG  
\$218.00  
**CONTOURED TORSO PAD**  
18.5 x 17 x 5



#YCEW  
\$265.00  
**PHILIPS CT HEADREST**  
12 x 12 x 6.5



#YCET  
\$229.00  
**CT/CATH LAB HEADREST**  
12.75 x 13 x 4.5  
(Neck Cradle 3" h)



#YCEV  
\$278.00  
**CONVERTIBLE CT HEADREST**  
14 x 13.5 x 6.75  
(2-Pc Flat and Convex Bases)


## Decubitus Block with Image Receptor Plate/Cassette Holder

Aides the technologist in positioning patients for cross-table lateral examinations. The decubitus pad ensures that the patient is elevated into a position above the plane of the imaging receptor for adequate imaging of the targeted body part. Provides excellent support and spine correction for laterals and decubitus views. Coated, firm foam.


#YCFQ | \$418.00

### Dimensions

| Lower Half | Back     | Channel   |
|------------|----------|-----------|
| 4" tall    | 11" tall | 0.5" wide |
| 36" long   | 36" long | 36" long  |
| 16" deep   | 10" deep | 3" deep   |



Convenient channel holds many types of image receptor plates and cassettes!



## Coated-Only Specialty Sponge Kits

### Specialty Coated CT Imaging Kit 4

#YSCT4 | \$763.00



### INCLUDED SPONGES


| QTY | INCLUDES | SPONGE TYPE      |
|-----|----------|------------------|
| (1) | YCAJ     | CT Knee Wedge    |
| (1) | YCBC     | 45° Wedge        |
| (1) | YCBF     | 27° Wedge        |
| (1) | YCBG     | 21° Wedge        |
| (1) | YCBH     | 16° Wedge Set    |
| (1) | YCBJ     | 12° Wedge Set    |
| (1) | YCFD     | 45° Spinal Wedge |

**Stocked sponges are marked Blue!**

## Specialty Sponge Kits

|  | COATED                   | NON-COATED               | INCLUDED SPONGES  |  |   |
|--|--------------------------|--------------------------|---|--|---|
|   | <b>#YSCA</b><br>\$263.00 | <b>#YKCA</b><br>\$63.00  | <b>coated</b><br>YCBB<br>YCAA<br>YCAK<br>YCBI                                 | <b>CLINIC KIT A</b><br><b>description</b><br>45° Wedge<br>Adult Obl Finger Block<br>Circular Head Block<br>15° Wedge Set   | <b>uncoated</b><br>YFBB<br>YFAA<br>YFAK<br>YFBI                                 |
|   | <b>#YSPA</b><br>\$403.00 | <b>#YKPA</b><br>\$155.00 | <b>coated</b><br>YCAC<br>YCAE<br>YCAK<br>YCAN<br>YCAQ<br>YCAX<br>YCBT         | <b>PEDIATRIC KIT</b><br><b>description</b><br>Bilateral Child Obl Finger Block<br>Pediatric Lateral Finger Block<br>Circular Head Block<br>AP/PA Child Head Holder<br>Child Head Immobilizer<br>Child Myelogram Block<br>Rectangle | <b>uncoated</b><br>YFAC<br>YFAE<br>YFAK<br>YFAN<br>YFAQ<br>YFAX<br>YFBT         |
|  | <b>#YSCB</b><br>\$544.00 | <b>#YKCB</b><br>\$132.00 | <b>coated</b><br>YCAA<br>YCAK<br>YCBC<br>YCBE<br>YCBI<br>YCBJ<br>YCBL<br>YCBD | <b>CLINIC KIT B</b><br><b>description</b><br>Adult Obl Finger Block<br>Circular Head Block<br>45° Wedge<br>34° Wedge Set<br>15° Wedge Set<br>12° Wedge Set<br>Rectangle<br>45° Wedge   | <b>uncoated</b><br>YFAA<br>YFAK<br>YFBC<br>YFBE<br>YFBI<br>YFBJ<br>YFBL<br>YFBD |

## General Sponge Kits

|   | COATED                   | NON-COATED              | INCLUDED SPONGES  |  |   |
|---|--------------------------|-------------------------|---|--|---|
|  | <b>#YSGP</b><br>\$297.00 | <b>#YKGP</b><br>\$77.00 | <b>coated</b><br>YCBD<br>YCBI<br>YCAK<br>YCAA                         | <b>GENERAL KIT P</b><br><b>description</b><br>45° Spine Wedge<br>15° Wedge<br>Circular Head Block<br>Oblique Finger Block  | <b>uncoated</b><br>YFBD<br>YFBI<br>YFAK<br>YFAA                         |
|  | <b>#YSGA</b><br>\$300.00 | <b>#YKGA</b><br>\$70.00 | <b>coated</b><br>YCAK<br>YCBC<br>YCBE<br>YCBI<br>YCBL                 | <b>GENERAL KIT A</b><br><b>description</b><br>Circular Head Block<br>45° Wedge<br>34° Wedge Set<br>15° Wedge Set<br>Rectangle  | <b>uncoated</b><br>YFAK<br>YFBC<br>YFBE<br>YFBI<br>YFBL                 |
|  | <b>#YSGB</b><br>\$405.00 | <b>#YKGB</b><br>\$96.00 | <b>coated</b><br>YCAA<br>YCAK<br>YCBC<br>YCBE<br>YCBI<br>YCBJ<br>YCBL | <b>GENERAL KIT B</b><br><b>description</b><br>Adult Obl Finger Block<br>Circular Head Block<br>45° Wedge<br>34° Wedge Set<br>15° Wedge Set<br>12° Wedge Set<br>Rectangle | <b>uncoated</b><br>YFAA<br>YFAK<br>YFBC<br>YFBE<br>YFBI<br>YFBJ<br>YFBL |


## Coated-Only Specialty Sponge Kits



COATED

INCLUDED SPONGES

### KNEE ARTHROPLASTY STRESS ASSEMBLY KIT



#YSAS  
\$334.00

45° Wedge (9x24x9)  
Contoured Cube (11x9x6)  
Non-Stealth


### ULTRASOUND SPECIALTY BREAST OR PARACENTESIS



#YSU8  
\$382.00

YCBG - 21° Wedge  
YCFD - 45° Spinal Wedge  
YCBU - 30°-60°-90° Multi-Angle Wedge

### ULTRASONOGRAPHER SCANNING KIT 3



#YSRC  
\$280.00

YCBT - (10) Individual Non-Stealth Rectangles

### GENERAL / PERINEAL KIT




#YSRD  
\$572.00

YCAV - Angular Bolster  
YCEE - 45° Echocardiography Wedge  
YCUZ - Echo-Ultrasound Wedge

## General Sponge Kits





|   | COATED                   | NON-COATED               | INCLUDED SPONGES   |  |  |
|---|--------------------------|--------------------------|--|--|--|
|    | <b>#YSGP</b><br>\$297.00 | <b>#YKGP</b><br>\$77.00  | <b>coated</b>  | <b>GENERAL KIT P</b>   | <b>uncoated</b>  |
|   |                          |                          | YCBD<br>YCBI<br>YCAK<br>YCAA   | <b>description</b><br>45° Spine Wedge<br>15° Wedge<br>Circular Head Block<br>Oblique Finger Block  | YFBD<br>YFBI<br>YFAK<br>YFAA   |
|    | <b>#YSGA</b><br>\$300.00 | <b>#YKGA</b><br>\$70.00  | <b>coated</b>  | <b>GENERAL KIT A</b>   | <b>uncoated</b>  |
|   |                          |                          | YCAK<br>YCBC<br>YCBE<br>YCBI<br>YCBL                                 | <b>description</b><br>Circular Head Block<br>45° Wedge<br>(2) 34° Wedges<br>(2) 15° Wedges<br>Rectangle  | YFAK<br>YFBC<br>YFBE<br>YFBI<br>YFBL                                 |
|   | <b>#YSGB</b><br>\$405.00 | <b>#YKGB</b><br>\$96.00  | <b>coated</b>  | <b>GENERAL KIT B</b>   | <b>uncoated</b>  |
|   |                          |                          | YCAA<br>YCAK<br>YCBC<br>YCBE<br>YCBI<br>YCBJ<br>YCBL                 | <b>description</b><br>Adult Obl Finger Block<br>Circular Head Block<br>45° Wedge<br>34° Wedge Set<br>15° Wedge Set<br>12° Wedge Set<br>Rectangle   | YFAA<br>YFAK<br>YFBC<br>YFBE<br>YFBI<br>YFBJ<br>YFBL                 |
|  | <b>#YSGD</b><br>\$438.00 | <b>#YKGD</b><br>\$126.00 | <b>coated</b>  | <b>GENERAL KIT D</b>   | <b>uncoated</b>  |
|   |                          |                          | YCAL<br>YCAN<br>YCAR<br>YCBF<br>YCBJ<br>YCBL                         | <b>description</b><br>Head & Neck Support<br>AP/PA Child Head Holder<br>Cervical Head Rest<br>27° Wedge<br>12° Wedge Set<br>Rectangle  | YFAL<br>YFAN<br>YFAR<br>YFBF<br>YFBJ<br>YFBL                         |
|  | <b>#YSGX</b><br>\$532.00 | <b>#YKGX</b><br>\$236.00 | <b>coated</b>  | <b>GENERAL KIT X</b>   | <b>uncoated</b>  |
|   |                          |                          | YCAK<br>YCBC<br>YCBE<br>YCBJ<br>YCBL<br>YCBN<br>Pad-51               | <b>description</b><br>Circular Head Block<br>45° Wedge<br>34° Wedge Set<br>12° Wedge Set<br>Rectangle<br>Rectangle<br>Standard Table Pad   | YFAK<br>YFBC<br>YFBE<br>YFBJ<br>YFBL<br>YFBN<br>Pad-51               |
|  | <b>#YSGE</b><br>\$849.00 | <b>#YKGE</b><br>\$190.00 | <b>coated</b>  | <b>GENERAL KIT E</b>   | <b>uncoated</b>  |
|   |                          |                          | YCAA<br>YCAK<br>YCAM<br>YCBB<br>YCBG<br>YCBI<br>YCBL<br>YCBN<br>YCBU | <b>description</b><br>Adult Obl Finger Block<br>Circular Head Block<br>AP/PA Adult Head Holder<br>45° Wedge<br>21° Wedge<br>15° Wedge Set<br>Rectangle<br>Rectangle<br>Multi-Angle Wedge | YFAA<br>YFAK<br>YFAM<br>YFBB<br>YFBG<br>YFBI<br>YFBL<br>YFBN<br>YFBU |



### General Sponge Kits



|   | COATED                   | NON-COATED               | INCLUDED SPONGES  |   |   |
|---|--------------------------|--------------------------|---|---|---|
|    | <b>#YSRA</b><br>\$182.00 | <b>#YKRA</b><br>\$38.00  | <b>coated</b><br>YCBQ<br>YCBR<br>YCBS<br>YCBT                                 | <b>RECTANGLE KIT</b><br><b>description</b><br>(3) Rectangles<br>(3) Rectangles<br>(3) Rectangles<br>(3) Rectangles  | <b>uncoated</b><br>YFBQ<br>YFBR<br>YFBS<br>YFBT                                 |
|    | <b>#YSWA</b><br>\$317.00 | <b>#YKWA</b><br>\$72.00  | <b>coated</b><br>YCBC<br>YCBF<br>YCBG<br>YCBI                                 | <b>WEDGE KIT A</b><br><b>description</b><br>45° Wedge<br>27° Wedge<br>21° Wedge<br>(2) 15° Wedges   | <b>uncoated</b><br>YFBC<br>YFBF<br>YFBG<br>YFBI                                 |
|   | <b>#YSWB</b><br>\$320.00 | <b>#YKWB</b><br>\$87.00  | <b>coated</b><br>YCBD<br>YCBE<br>YCBF<br>YCBG                                 | <b>WEDGE KIT B</b><br><b>description</b><br>45° Spine Wedge<br>(2) 34° Wedges<br>27° Wedge<br>21° Wedge   | <b>uncoated</b><br>YFBD<br>YFBE<br>YFBF<br>YFBG                                 |
|  | <b>#YSOA</b><br>\$408.00 | <b>#YKOA</b><br>\$120.00 | <b>coated</b><br>YCAA<br>YCAM<br>YCAR<br>YCBD<br>YCBL                         | <b>ORTHO KIT A</b><br><b>description</b><br>Adult Obl Finger Block<br>AP/PA Adult Head Holder<br>Cervical Head Rest<br>45° Spine Wedge<br>Rectangle   | <b>uncoated</b><br>YFAA<br>YFAM<br>YFAR<br>YFBD<br>YFBL                         |
|  | <b>#YSOB</b><br>\$518.00 | <b>#YKOB</b><br>\$131.00 | <b>coated</b><br>YCAL<br>YCAN<br>YCAR<br>YCBC<br>YCBG<br>YCBJ<br>YCBL         | <b>ORTHO KIT B</b><br><b>description</b><br>Head & Neck Support<br>AP/PA Child Head Holder<br>Cervical Head Rest<br>45° Wedge<br>21° Wedge<br>(2) 12° Wedges<br>Rectangle   | <b>uncoated</b><br>YFAL<br>YFAN<br>YFAR<br>YFBC<br>YFBG<br>YFBJ<br>YFBL         |
|  | <b>#YSOC</b><br>\$794.00 | <b>#YKOC</b><br>\$226.00 | <b>coated</b><br>YCAA<br>YCAD<br>YCAF<br>YCAH<br>YCAK<br>YCAL<br>YCAR<br>YCBU | <b>ORTHO KIT C</b><br><b>description</b><br>Adult Obl Finger Block<br>Adult Lateral Finger Block<br>Leg Immobilizer<br>Knee Arthrogram Block<br>Circular Head Block<br>Head & Neck Support<br>Cervical Head Rest<br>Multi-Angle Wedge | <b>uncoated</b><br>YFAA<br>YFAD<br>YFAF<br>YFAH<br>YFAK<br>YFAL<br>YFAR<br>YFBU |

## Stocked Vinyl-Covered Positioning Bolsters

Reusable blocks, bolsters, and wedges are manufactured from firm, high-density urethane foam. Royal Blue outer covering is durable, heavy gauge, reinforced vinyl. Seams are sewn and tear resistant. Easy to clean Bolsters are **RADIOLUCENT**. Seams may show. Provides excellent support for a wide variety of imaging needs.

### STOCKED BOLSTERS

| CAT.# | Description             | Size (W x L x H)    | Price    |
|-------|-------------------------|---------------------|----------|
| YVF40 | 26° Knee Wedge          | 13" x 17" x 6"      | \$320.00 |
| YVFC  | 20° Wedge               | 24" x 18" x 6"      | \$193.00 |
| YVFD  | 45° Wedge               | 8" x 16" x 8"       | \$219.00 |
| YVFE  | 20° Incline Wedge       | 20" x 32" x 12"     | \$277.00 |
| YVFM  | Knee Wedge              | 16" x 28.2" x 8.25" | \$370.00 |
| YVFO  | 45° Spinal Wedge        | 7" x 21" x 7"       | \$248.00 |
| YVFU  | 30°/60°/90° Wedge       | 10" x 24" x 6"      | \$227.00 |
| YVJV  | Comfort Pillow w/ Strap | 16" x 25" x 8.25"   | \$288.00 |
| YVFZ  | Endo-Ultrasound Wedge   | 20" x 22" x 6"      | \$484.00 |
| YVFG  | Rectangle               | 8" x 12" x 16"      | \$106.00 |
| YVFH  | Rectangle               | 12" x 14" x 3"      | \$72.00  |
| YVFI  | Rectangle               | 24" x 15" x 3"      | \$127.00 |
| YVFJ  | Rectangle               | 8" x 12" x 3"       | \$81.00  |
| YVFK  | Rectangle               | 20" x 13" x 3"      | \$118.00 |
| YVFP  | Decubitus Block         | 18" x 24" x 4"      | \$227.00 |
| YVFQ  | DEXA Block              | 20" x 20" x 12"     | \$351.00 |
| YVFB  | Semi-Cylinder           | 8" x 18" x 6"       | \$127.00 |

**Stocked sponges are ready to ship!**



## Body Comfort Wedge Pillow

High density foam enclosed in a durable, heavy gauge, reinforced Royal Blue vinyl covering. Supports patient's head and shoulders when arms must be held above head. Dual hand grips reduce patient fatigue. Measures: 16" x 25" x 8.25"

#YVJV | \$288.00



**Easy access grips reduces patient fatigue**

## Non-Stocked Vinyl-Covered Positioning Bolsters

## NON-STOCKED WEDGE BOLSTERS

| CAT.#  | Description      | Size (W x L x H)       | Price    |
|--------|------------------|------------------------|----------|
| YVF155 | 13.5° Knee Wedge | 10" x 10" x 3"         | \$450.00 |
| YVBH   | 16° Knee Wedge   | 8.5" x 11" x 3.5"      | \$287.00 |
| YVBG   | 21° Knee Wedge   | 12.25" x 14.25" x 4.5" | \$329.00 |
| YVF156 | 21° Knee Wedge   | 10" x 10" x 4.6"       | \$455.00 |
| YVDB   | 30° Knee Wedge   | 28" x 27" x 16"        | \$366.00 |
| YVDC   | 30° Knee Wedge   | 9" x 30" x 9"          | \$221.00 |
| YVF126 | 30° Knee Wedge   | 23" x 11" x 6.5"       | \$196.00 |
| YVDA   | 45° Knee Wedge   | 14" x 28" x 14"        | \$220.00 |
| YVF129 | 45° Knee Wedge   | 5" x 15" x 5"          | \$103.00 |
| YVF178 | 45° Knee Wedge   | 9" x 30" x 9"          | \$436.00 |

## NON-STOCKED SQUARE BOLSTERS

| CAT.#  | Description | Size (W x L x H) | Price    |
|--------|-------------|------------------|----------|
| YVF172 | Square      | 18" x 18" x 4"   | \$369.00 |
| YVF42  | Square      | 12" x 12" x 9"   | \$539.00 |

## NON-STOCKED WEDGE BOLSTERS

| CAT.#  | Description   | Size (W x L x H)        | Price    |
|--------|---------------|-------------------------|----------|
| YVF161 | Cylinder      | 4.5" x 18" x 4.5"       | \$384.00 |
| YVF167 | Cylinder      | 4" x 48" x 4"           | \$152.00 |
| YVF53  | Cylinder      | 4.5" x 12" x 4.5"       | \$401.00 |
| YVFA   | Cylinder      | 8" x 24" x 8"           | \$155.00 |
| YVFS   | Cylinder      | 6" x 17" x 6"           | \$193.00 |
| YVF119 | Semi-Cylinder | 12.75" x 22.75" x 8.75" | \$293.00 |
| YVF164 | Semi-Cylinder | 24" x 36" x 16"         | \$695.00 |
| YVF180 | Semi-Cylinder | 6" x 12" x 4.5"         | \$418.00 |
| YVF46  | Semi-Cylinder | 14" x 20" x 5"          | \$375.00 |

**Stocked sponges  
are marked Blue!**

## NON-STOCKED RECTANGLE BOLSTERS

| CAT.#  | Description | Size (W x L x H)  | Price    |
|--------|-------------|-------------------|----------|
| YVF110 | Rectangle   | 6" x 12" x 1"     | \$105.00 |
| YVF109 | Rectangle   | 6.5" x 9" x 1.5"  | \$118.00 |
| YVF137 | Rectangle   | 5" x 11" x 1.5"   | \$142.00 |
| YVF133 | Rectangle   | 4.25" x 34" x 2"  | \$144.00 |
| YVF160 | Rectangle   | 6" x 8" x 2"      | \$316.00 |
| YVF73  | Rectangle   | 14" x 24" x 2"    | \$112.00 |
| YVF75  | Rectangle   | 4.5" x 22" x 2"   | \$87.00  |
| YVF78  | Rectangle   | 17" x 25" x 2"    | \$420.00 |
| YVF79  | Rectangle   | 7.5" x 8.5" x 2"  | \$127.00 |
| YVF93  | Rectangle   | 24" x 36" x 2"    | \$278.00 |
| YVF170 | Rectangle   | 12" x 15" x 2.5"  | \$338.00 |
| YVF60  | Rectangle   | 14" x 23" x 2.5"  | \$138.00 |
| YVF74  | Rectangle   | 13" x 22" x 2.5"  | \$364.00 |
| YVF105 | Rectangle   | 14" x 18" x 3"    | \$350.00 |
| YVF136 | Rectangle   | 5" x 11" x 3"     | \$214.00 |
| YVF177 | Rectangle   | 12" x 4" x 3"     | \$326.00 |
| YVF107 | Rectangle   | 20" x 24" x 4"    | \$280.00 |
| YVF108 | Rectangle   | 20" x 7.5" x 4"   | \$180.00 |
| YVF120 | Rectangle   | 18" x 36" x 4"    | \$289.00 |
| YVF134 | Rectangle   | 5" x 30" x 4"     | \$196.00 |
| YVF57  | Rectangle   | 3" x 14" x 17"    | \$145.00 |
| YVF179 | Rectangle   | 4" x 10" x 4"     | \$315.00 |
| YVF182 | Rectangle   | 9.5" x 9.5" x 4"  | \$372.00 |
| YVF184 | Rectangle   | 10" x 20" x 4"    | \$386.00 |
| YVF102 | Rectangle   | 8" x 10" x 6"     | \$155.00 |
| YVF132 | Rectangle   | 12" x 118" x 6"   | \$249.00 |
| YVF175 | Rectangle   | 14" x 14" x 6"    | \$403.00 |
| YVF186 | Rectangle   | 8" x 20" x 6"     | \$385.00 |
| YVF188 | Rectangle   | 10" x 20" x 6"    | \$396.00 |
| YVF169 | Rectangle   | 6.5" x 21" x 6.5" | \$341.00 |
| YVFT   | Rectangle   | 16" x 12" x 7"    | \$331.00 |
| YVF77  | Rectangle   | 6" x 18" x 7.5"   | \$320.00 |
| YVF56  | Rectangle   | 2" x 8" x 10"     | \$106.00 |
| YVF99  | Rectangle   | 14" x 17" x 10"   | \$476.00 |
| YVF187 | Rectangle   | 16" x 16" x 12"   | \$463.00 |
| YVF43  | Rectangle   | 12" x 12" x 12"   | \$611.00 |
| YVF84  | Rectangle   | 14" x 17" x 12"   | \$446.00 |
| YVFW   | Rectangle   | 16" x 18" x 12"   | \$383.00 |
| YVF176 | Rectangle   | 14.5" x 18" x 13" | \$480.00 |

## Vinyl-Covered Positioning Bolster Kits

Reusable blocks, bolsters, and wedges are manufactured from firm, high-density urethane foam. Royal Blue outer covering is durable, heavy gauge, reinforced vinyl. Seams are sewn and tear resistant. Easy to clean Bolsters are RADIOLUCENT. Seams may show. Provides excellent support for a wide variety of imaging needs.

### GENERAL KIT 1 | YSU1 - \$399.00

| CAT.# | Description   | Size (W x L x H) | QTY |
|-------|---------------|------------------|-----|
| YVFB  | Semi-Cylinder | 8" x 18" x 6"    | 1   |
| YVFG  | Rectangle     | 8" x 12" x 16"   | 2   |
| YVFJ  | Rectangle     | 8" x 12" x 3"    | 1   |

### GENERAL KIT 2 | YSU2 - \$275.00

| CAT.# | Description | Size (W x L x H) | QTY |
|-------|-------------|------------------|-----|
| YVFG  | Rectangle   | 8" x 12" x 16"   | 2   |
| YVFJ  | Rectangle   | 8" x 12" x 3"    | 1   |

### GENERAL KIT 3 | YSU3 - \$429.00

| CAT.# | Description   | Size (W x L x H) | QTY |
|-------|---------------|------------------|-----|
| YVFB  | Semi-Cylinder | 8" x 18" x 6"    | 1   |
| YVFD  | 45° Wedge     | 8" x 16" x 8"    | 1   |
| YVFG  | Rectangle     | 8" x 12" x 16"   | 1   |

### GENERAL KIT 4 | YSU4 - \$1015.00

| CAT.# | Description   | Size (W x L x H) | QTY |
|-------|---------------|------------------|-----|
| YVFA  | Cylinder      | 8" x 24" x 8"    | 1   |
| YVFB  | Semi-Cylinder | 8" x 18" x 6"    | 1   |
| YVFC  | 20° Wedge     | 24" x 18" x 6"   | 1   |
| YVFD  | 45° Wedge     | 8" x 16" x 8"    | 1   |
| YVFG  | Rectangle     | 8" x 12" x 16"   | 1   |
| YVFH  | Rectangle     | 12" x 14" x 3"   | 1   |
| YVFJ  | Rectangle     | 8" x 12" x 3"    | 1   |
| YVFK  | Rectangle     | 20" x 13" x 3"   | 1   |

### GENERAL KIT 5 | YSU5 - \$825.00

| CAT.# | Description      | Size (W x L x H) | QTY |
|-------|------------------|------------------|-----|
| YVFB  | Semi-Cylinder    | 8" x 18" x 6"    | 1   |
| YVFC  | 20° Wedge        | 24" x 18" x 6"   | 1   |
| YVFD  | 45° Wedge        | 8" x 16" x 8"    | 1   |
| YVFJ  | Rectangle        | 8" x 12" x 3"    | 1   |
| YVFO  | 45° Spinal Wedge | 7" x 21" x 7"    | 1   |

### GENERAL KIT 6 | YSU6 - \$775.00

| CAT.# | Description      | Size (W x L x H) | QTY |
|-------|------------------|------------------|-----|
| YVFJ  | Rectangle        | 8" x 12" x 3"    | 1   |
| YVFK  | Rectangle        | 20" x 13" x 3"   | 1   |
| YVFO  | 45° Spinal Wedge | 7" x 21" x 7"    | 1   |
| YVFM  | Knee Wedge       | 16" x 28.2" x 8" | 1   |

**Don't see the kit you're looking for? We are happy to help build your perfect solution!**  
**Contact us today**



## Vinyl-Covered Specialty Bolsters

### KNEE WEDGE

#YVFM | \$370.00

16" x 25" x 8.25"



### ENDO ULTRASOUND VINYL WEDGE

#YVFZ | \$484.00

20" x 23" x 6"



### WHEELCHAIR PATIENT POSITIONING BOLSTER

This patient positioning bolster is perfect for propping up cassettes behind patients in wheelchairs. Sponge is curved to fit in most wheelchairs. Covered with heavy duty vinyl royal blue vinyl.

#YVWS | \$595.00

18" x 20" x 5"



### TAPERED STRETCHER/WHEELCHAIR BOLSTER

This patient positioning bolster is perfect for stretcher and wheelchair imaging and mobile x-ray. Covered with heavy duty vinyl royal blue vinyl.

#YVWA | \$759.00

14"(W) x 24.5"(L) x 1.5"-5.5"(H)



## Vinyl-Covered Specialty Bolster Kits

### DEXA POSITIONING KIT

The large bolster aids in positioning the AP lumbar spine. The small head bolster and 15° "A" frame splint lock the feet into position to abduct both legs to hold full extension at the knees and hips.

#TI-SET | \$608.00

| CAT# | Description      | Size (W x L x H) | QTY |
|------|------------------|------------------|-----|
| YVFJ | Head Bolster     | 8" x 12" x 3"    | 1   |
| YVFQ | DEXA Leg Bolster | 20" x 20" x 12"  | 1   |

| Description                       | Size (W x T x D) | QTY |
|-----------------------------------|------------------|-----|
| "A" Frame 15° Locking Foot Splint | 7.5" x 12" x 8"  | 1   |



### LATERAL DECUBITUS POSITIONER KIT

This KIT includes a 90° support with hook and loop stability straps and a YVFK decubitus bolster to aid in cross lateral examinations. The pad assures the patient is elevated into the imaging field.

#TBC-2214 | \$588.00

Pad Size: 20"(W) x 13"(L) x 3"(H)







ADMINISTRATION/MANUFACTURING

7117 Centennial Blvd. Nashville, TN 37209

Orders can be placed online or sent to [orders@techno-aide.com](mailto:orders@techno-aide.com)

For official Quotes: [quotes@techno-aide.com](mailto:quotes@techno-aide.com)

Any questions can be directed to [customerservice@techno-aide.com](mailto:customerservice@techno-aide.com)

P: 800.251.2629 or 615.350.7030

F: 800.554.6275 or 615.350.7879

[techno-aide.com](http://techno-aide.com)