



# Table Pads

2020  
PRODUCT  
GUIDE

## Picking the Perfect Table Pad

Techno-Aide is here to help you build your perfect table pad! Our table pads are specially designed to meet a variety of imaging needs for nearly every x-ray table! Below is a simplified process for table pad ordering.\*

### TABLE PAD FOAM TYPES

#### COMFORT FOAM

**Our most popular foam!** The perfect foam for patient comfort and x-ray imaging. Available in different sizes and thicknesses.



1.8lb | 60 ILD

#### FIRM DENSITY FOAM

**Perfect for longer procedures!** This high density foam provides more patient support than our comfort foam. Perfect for extended imaging.



2.0lb | 65 ILD

#### MEMORY FOAM

Great for long procedures and sensitive patients. 1" of 5.5lb memory foam is placed on top of 1" of comfort foam.



5.5lb | 65 ILD

### TABLE PAD COVERINGS

#### UPHOLSTERY VINYL

**Our most popular cover!** Comfortable, durable, heavy gauge, medical-grade, reinforced vinyl. [1 mil.]

See PG. 11 for more color options!



Black Upholstery  
Default Color



Green Upholstery



Blue Upholstery

Features sewn tapered ends


#### NYLON FABRIC

A lightweight, filament reinforced nylon fabric.

Now available in feather or diamond ripstop! PG. 11 for color options!



Diamond Ripstop  
61-Sapphire  
Default Color



Feather Ripstop  
R710-Gray

#### ECONOMY VINYL

A lightweight, inexpensive, light blue vinyl.

Add ons/upgrades not available.  
Will have fused tapered ends.



Light Blue  
Economy Vinyl

#### TECHNO-TUFF™

Heavyweight, mesh, laminate, reinforced vinyl.

Features handles for easy patient transport!



TT Black  
Default Color



TT Blue

#### 4-WAY STRETCH

A comfortable, pliable, anti-microbial premium fabric.

Only available with premium specialty pads.




Black 4-Way Stretch

### TABLE PAD END FINISHES

#### TAPERED ENDS

Standard closure with sewn or fused side seams.

Includes grommets for proper table pad storage.



#### BOXED ENDS

Completely sewn covering with squared finished corners.

Only available with premium specialty pads!



\* Seams will show on all table pads under normal imaging.

## Radiolucent X-Ray Table Pads

Pads can be used for imaging tables or stretchers. Manufactured for durability and longevity. Helps maintain complete control of the patient. Affordable for each room and table. Standard table pads come with tapered ends and two brass eyelets at each end for hanging. Cleans easy with mild cleaner or germicide.



### STEP 1 - CHOOSE A FOAM TYPE

Comfort Foam - Our most popular table pad foam.

Firm Foam - Higher density for more patient support.

### STEP 2 - CHOOSE THE BEST FIT

#### STANDARD

##### Comfort Foam - 72" x 23.25" x 1"

CAT.#	Color	Fabric	Price
PAD-51	Light Blue	Econo Vinyl	\$118.00
PAD-53	Green	Upholstery Vinyl	\$208.00
PAD-54	Red	Upholstery Vinyl	\$208.00
PAD-55	Black	Upholstery Vinyl	\$208.00
PAD-50	Blue	Upholstery Vinyl	\$208.00
PAD-52	Yellow	Upholstery Vinyl	\$208.00
PAD-80	See Page 11	Nylon Fabric	\$178.00

##### Firm Foam - 72" x 23.25" x 2"

CAT.#	Color	Fabric	Price
PAH-51	Light Blue	Standard Vinyl	\$178.00
PAH-53	Green	Upholstery Vinyl	\$268.00
PAH-54	Red	Upholstery Vinyl	\$268.00
PAH-55	Black	Upholstery Vinyl	\$268.00
PAH-50	Blue	Upholstery Vinyl	\$268.00
PAH-52	Yellow	Upholstery Vinyl	\$268.00
PAH-80	See Page 11	Nylon Fabric	\$228.00

#### STANDARD PLUS

##### Comfort Foam - 80" x 30" x 1"

CAT.#	Color	Fabric	Price
PAD-SP51	Light Blue	Standard Vinyl	\$228.00
PAD-SP53	Green	Upholstery Vinyl	\$318.00
PAD-SP54	Red	Upholstery Vinyl	\$318.00
PAD-SP55	Black	Upholstery Vinyl	\$318.00
PAD-SP50	Blue	Upholstery Vinyl	\$318.00
PAD-SP52	Yellow	Upholstery Vinyl	\$318.00
PAD-SP80	See Page 11	Nylon Fabric	\$288.00

##### Firm Foam - 80" x 30" x 2"

CAT.#	Color	Fabric	Price
PAH-SP51	Light Blue	Standard Vinyl	\$288.00
PAH-SP53	Green	Upholstery Vinyl	\$368.00
PAH-SP54	Red	Upholstery Vinyl	\$368.00
PAH-SP55	Black	Upholstery Vinyl	\$368.00
PAH-SP50	Blue	Upholstery Vinyl	\$368.00
PAH-SP52	Yellow	Upholstery Vinyl	\$368.00
PAH-SP80	See Page 11	Nylon Fabric	\$338.00

#### STANDARD LONG

##### 77" x 23.25" x [Foam Type]

CAT.#	Foam Type	Thickness	Price
PAD-LS77	Comfort	1"	\$318.00
PAH-LS77	Firm	2"	\$378.00

#### STANDARD NARROW

##### 72" x 19" x [Foam Type]

CAT.#	Foam Type	Thickness	Price
PAD-SN19	Comfort	1"	\$289.00
PAH-SN19	Firm	2"	\$378.00

#### STANDARD LONG & NARROW

##### 77" x 19" x [Foam Type]

CAT.#	Foam Type	Thickness	Price
PAD-LN79	Comfort	1"	\$318.00
PAH-LN79	Firm	2"	\$378.00


### STEP 3 - CHOOSE A PAD COVERING

#### UPHOLSTERY VINYL

Durable, Heavy Gauge, Medical-Grade Vinyl, Reinforced 1 mil.




Black  
Upholstery



Green  
Upholstery



Red  
Upholstery



Blue  
Upholstery



Yellow  
Upholstery

#### ECONO-VINYL

Lightweight, Economical Vinyl



Light Blue  
(Economy vinyl)

\* Seams will show on all table pads under normal imaging.

## Radiolucent Premium Specialty Pads

These premium table pads offer excellent support and durability. All premium specialty pads feature squared off corners that provide an exceptional finished look. Available in blue or black upholstery vinyl and our premium 4-way stretch health-care fabric. Pads clean easily with a mild cleaner or germicide.



Fully sewn  
4-way cover  
with squared  
corners!

### PREMIUM PAD

#### Comfort Foam - 72" x 23.25" x 1"

CAT.#	Color	Covering	Price
PPAD-75	Black	4-Way Stretch	\$498.00
PPAD-55	Black	Upholstery Vinyl	\$438.00
PPAD-50	Blue	Upholstery Vinyl	\$438.00

#### Firm Foam - 72" x 23.25" x 2"

CAT.#	Color	Covering	Price
PPAH-75	Black	4-Way Stretch	\$678.00
PPAH-55	Black	Upholstery Vinyl	\$468.00
PPAH-50	Blue	Upholstery Vinyl	\$468.00

#### Memory Foam\* - 72" x 23.25" x 2"

CAT.#	Color	Covering	Price
PPAM-80	Black	4-Way Stretch	\$658.00

\*these pads are 1" of bottom layer standard foam and 1" of top layer 5.5lb memory foam

### PREMIUM PLUS PAD

#### Comfort Foam - 80" x 30" x 1"

CAT.#	Color	Covering	Price
PPAD-SP75	Black	4-Way Stretch	\$498.00
PPAD-SP55	Black	Upholstery Vinyl	\$438.00
PPAD-SP50	Blue	Upholstery Vinyl	\$438.00

#### Firm Foam - 80" x 30" x 2"

CAT.#	Color	Covering	Price
PPAH-SP75	Black	4-Way Stretch	\$498.00
PPAH-SP55	Black	Upholstery Vinyl	\$438.00
PPAH-SP50	Blue	Upholstery Vinyl	\$438.00

#### Memory Foam\* - 80" x 30" x 2"

CAT.#	Color	Covering	Price
PPAM-SP80	Black	4-Way Stretch	Call for Price

\*these pads are 1" of bottom layer standard foam and 1" of top layer 5.5lb memory foam

## Radiolucent Mobility Table Pads

Our mobility table pads feature Techno-Tuff™ vinyl. Techno-Tuff is a low-friction, waterproof, stain proof, antimicrobial, durable mesh, laminate reinforced vinyl. This material includes heavy-duty handles along both sides for moving patients from a stretcher to a table. Not intended to lift patients. Manufactured to transfer patients laterally from surface to surface.


### STANDARD MOBILITY PAD

#### Comfort Foam - 72" x 23.25" x 1"

CAT.#	Color	Fabric	Price
PAD-H8R	Royal Blue	Techno-Tuff™	\$298.00
PAD-H8K	Black	Techno-Tuff™	\$298.00

#### Firm Foam - 72" x 23.25" x 2"

CAT.#	Color	Fabric	Price
PAH-H8R	Royal Blue	Techno-Tuff™	\$358.00
PAH-H8K	Black	Techno-Tuff™	\$358.00



### STANDARD PLUS MOBILITY PAD

#### Comfort Foam - 80" x 30" x 1"

CAT.#	Color	Fabric	Price
PAD-SPH8R	Royal Blue	Techno-Tuff™	\$375.00
PAD-SPH8K	Black	Techno-Tuff™	\$375.00

#### Firm Foam - 80" x 30" x 2"

CAT.#	Color	Fabric	Price
PAH-SPH8R	Royal Blue	Techno-Tuff™	\$425.00
PAH-SPH8K	Black	Techno-Tuff™	\$425.00

### TECHNO-TUFF VINYL

Heavyweight, Mesh Laminate,  
Anti-Microbial Reinforced Vinyl



TT Black



TT Blue

Extreme  
Durability

\* Seams will show on all table pads under normal imaging.

## Radiolucent Stretcher Replacement Pads

Economically priced replacement stretcher pads offer excellent support and durability. The black 4-way-stretch antimicrobial fabric is a tailored fit over our firm polyfoam. Squared off corners provide a finished look with a resealable hook and loop end panel. Pads are CAL 117 compliant and clean easily with mild cleaner or germicide.

All stretcher replacement pads are 3" in height.

### WITH ANTI-SLIP BACKING

#### 72" LENGTH

CAT. #	Length x Width	Price
PMR-7224ASB	72" x 24"	\$768.00
PMR-7226ASB	72" x 26"	\$788.00
PMR-7228ASB	72" x 28"	\$838.00
PMR-7230ASB	72" x 30"	\$888.00

#### 74" LENGTH

CAT. #	Length x Width	Price
PMR-7424ASB	74" x 24"	\$778.00
PMR-7426ASB	74" x 26"	\$798.00
PMR-7428ASB	74" x 28"	\$858.00
PMR-7430ASB	74" x 30"	\$948.00

#### 76" LENGTH

CAT. #	Length x Width	Price
PMR-7624ASB	76" x 24"	\$778.00
PMR-7626ASB	76" x 26"	\$808.00
PMR-7628ASB	76" x 28"	\$858.00
PMR-7630ASB	76" x 30"	\$918.00

### WITHOUT BACKING

#### 72" LENGTH

CAT. #	Length x Width	Price
PMR-7224	72" x 24"	\$638.00
PMR-7226	72" x 26"	\$648.00
PMR-7228	72" x 28"	\$698.00
PMR-7230	72" x 30"	\$748.00

#### 74" LENGTH

CAT. #	Length x Width	Price
PMR-7424	74" x 24"	\$648.00
PMR-7426	74" x 26"	\$658.00
PMR-7428	74" x 28"	\$718.00
PMR-7430	74" x 30"	\$768.00

#### 76" LENGTH

CAT. #	Length x Width	Price
PMR-7624	76" x 24"	\$648.00
PMR-7626	76" x 26"	\$668.00
PMR-7628	76" x 28"	\$718.00
PMR-7630	76" x 30"	\$778.00

### CORNER OPTIONS FOR STRETCHER REPLACEMENT PADS

#### 4 SQUARED CORNERS

##### #PMR

Default replacement stretcher pad with 4 squared corners. For stretchers without hand holes at any corner.



#### 2 SLANTED CORNERS

##### #PMR-A

Special replacement stretcher pad with 2 slanted corners. For stretchers with 2 hand holes at only one end.



#### 4 SLANTED CORNERS

##### #PMR-B

Special replacement stretcher pad with 4 slanted corners. For stretchers with 4 hand holes at both ends.




Custom sizes, including thickness, are available upon request. For more information, please email [CustomerService@Techno-Aide.com](mailto:CustomerService@Techno-Aide.com)

## Manufacturer Specific Table Pads

Techno-Aide is proud to present our extensive line of Manufacturer Specific table pads. These pads are designed to fit on specific manufacturers' imaging tables. Each pad offers excellent support and durability. Our premium, black, 4-way-stretch, antimicrobial fabric is tailor-fitted over your choice of either 2" premium comfort foam or our 2" two-ply foam with memory. Squared off corners and ends provide a nice, finished look that is perfect for any imaging table. Pads clean easily with mild cleaner or germicide.



Now available in both  
Comfort Foam and  
Memory Foam options!




Spectrum Angio Pad

Overall Dimensions  
90" x 30"

**COMFORT**  
#PPMS-ANS  
\$725.00

**MEMORY**  
#PMMS-ANS  
\$1115.00



Cysto Table Pad

Overall Dimensions  
44" x 24"

**COMFORT**  
#PPMS-CYS  
\$350.00

**MEMORY**  
#PMMS-CYS  
\$503.00



Memco X-Ray

Overall Dimensions  
83" x 19"

**COMFORT**  
#PPMS-MEM  
\$525.00

**MEMORY**  
#PMMS-MEM  
\$770.00



Shimadzu KS50

Overall Dimensions  
97" x 18"

**COMFORT**  
#PPMS-S50  
\$996.00

**MEMORY**  
#PMMS-S50  
\$665.00



XRE 8844A0

Overall Dimensions  
80" x 18.5"

**COMFORT**  
#PPMS-XE8  
\$625.00

**MEMORY**  
#PMMS-XE8  
\$860.00




XRE-ADC

Overall Dimensions  
94.5" x 28"


**COMFORT**  
#PPMS-XEA  
\$835.00

**MEMORY**  
#PMMS-XEA  
\$1335.00



**COMFORT**  
#PPGE-98R  
\$645.00

**MEMORY**  
#PMGE-98R  
\$955.00



**COMFORT**  
#PPGE-AFM  
\$960.00

**MEMORY**  
#PMGE-AFM  
\$1629.00



**COMFORT**  
#PPGE-98N  
\$525.00

**MEMORY**  
#PMGE-98N  
\$780.00



**COMFORT**  
#PPGE-ANG  
\$825.00

**MEMORY**  
#PMGE-ANG  
\$1326.00




**COMFORT**  
#PPGE-NRW  
\$630.00

**MEMORY**  
#PMGE-NRW  
\$842.00




**COMFORT**  
#PPGE-LUP  
\$950.00

**MEMORY**  
#PMGE-LUP  
\$1517.00



**COMFORT**  
#PPGE-OM3  
\$1045.00

**MEMORY**  
#PMGE-OM3  
\$1317.00



**COMFORT**  
#PPGE-OM4  
\$1264.00

**MEMORY**  
#PMGE-OM4  
\$1629.00



**COMFORT**  
#PPGE-OM5  
\$1100.00

**MEMORY**  
#PMGE-OM5  
\$1746.00

Now available in both  
Comfort Foam and  
Memory Foam options!

### Siemens Angio Pad

Overall Dimensions  
66" x 21"

**COMFORT**  
#PPSI-ANG  
\$630.00

**MEMORY**  
#PMSI-ANG  
\$862.00

### Siemens Koordinat Angio Table Pad

Overall Dimensions  
79" x 21"

**COMFORT**  
#PPSI-ANK  
\$690.00

**MEMORY**  
#PMSI-ANK  
\$1056.00

### Siemens CCL Cardiac Cath

Overall Dimensions  
79" x 23.5"

**COMFORT**  
#PPSI-CCL  
\$720.00

**MEMORY**  
#PMSI-CCL  
\$1098.00

### Siemens Orbiter ECT Special Procedure

Overall Dimensions  
78" x 18.5"

**COMFORT**  
#PPSI-ECT  
\$660.00

**MEMORY**  
#PMSI-ECT  
\$976.00

### Siemens HiQ

Overall Dimensions  
96" x 21"

**COMFORT**  
#PPSI-HIQ  
\$680.00

**MEMORY**  
#PMSI-HIQ  
\$1023.00

### Siemens Somatom DR

Overall Dimensions  
82" x 16"

**COMFORT**  
#PPSI-SDR  
\$520.00

**MEMORY**  
#PMSI-SDR  
\$746.00

### Siemens HiQ Short

Overall Dimensions  
76" x 15"

**COMFORT**  
#PPSI-HIS  
\$480.00

**MEMORY**  
#PMSI-HIS  
\$700.00

### Siemens Somatom Plus

Overall Dimensions  
75" x 16"

**COMFORT**  
#PPSI-SPL  
\$480.00

**MEMORY**  
#PMSI-SPL  
\$700.00

### Siemens HiQ XL

Overall Dimensions  
114" x 15"

**COMFORT**  
#PPSI-HIX  
\$640.00

**MEMORY**  
#PMSI-HIX  
\$888.00

Now available in both  
Comfort Foam and  
Memory Foam options!

### Toshiba Cat 150

Overall Dimensions  
91" x 24"

#### COMFORT

#PPTO-150

\$730.00

#### MEMORY

#PMTO-150

\$1120.00

### Toshiba Angio Table Pad

Overall Dimensions  
104" x 21"

#### COMFORT

#PPTO-ANG

\$790.00

#### MEMORY

#PMTO-ANG

\$1116.00

### Toshiba Angio Contoured

Overall Dimensions  
96" x 24"

#### COMFORT

#PPTO-ANC

\$730.00

#### MEMORY

#PMTO-ANC

\$1121.00

### Toshiba Cat 250

Overall Dimensions  
104" x 27.5"

#### COMFORT

#PPTO-250

\$985.00

#### MEMORY

#PMTO-250

\$1629.00

### Toshiba Special Procedure

Overall Dimensions  
94" x 21"

#### COMFORT

#PPTO-SCP

\$710.00

#### MEMORY

#PMTO-SCP

\$1079.00

### Toshiba Cat FX-2

Overall Dimensions  
112" x 26"

#### COMFORT

#PPTO-FX2

\$960.00

#### MEMORY

#PMTO-FX2

\$1638.00

### Picker Angio

Overall Dimensions  
99" x 22"

#### COMFORT

#PPPI-ANG

\$730.00

#### MEMORY

#PMPI-ANG

\$1116.00

### Picker 1200

Overall Dimensions  
82" x 16"

#### COMFORT

#PPPI-120

\$520.00

#### MEMORY

#PMPI-120

\$746.00

### Picker CT

Overall Dimensions  
89" x 16"

#### COMFORT

#PPPI-PCT

\$530.00

#### MEMORY

#PMPI-PCT

\$785.00

### Nuclear X-Ray

Overall Dimensions  
71" x 24"

#### COMFORT

#PPPI-NXR

\$640.00

#### MEMORY

#PMPI-NXR

\$916.00

## Table Pad Add-Ons

Techno-Aide offers multiple premium add-ons to help build your perfect table pad. Below are our most popular premium upgrades to choose from! **Table Pad add-ons are not available on PAD-51 Economy vinyl.**

### CHOOSE A STABILITY OPTION

#### ANTI-SKID BACKING

**#PAD-ASB | \$75.00** per pad

This add-on minimizes movement of table pads on any imaging table, but allows easy removal of your table pad.

**Not available on PAD-51 economy vinyl.**



Our most popular option for table pad stability!

OR

#### HOOK & LOOP

**#PAD-HKLP | \$45.00** per pad

This add-on securely places a table pad on any imaging table. Allows for removal of pad from table.

**Not available on PAD-51 economy vinyl.**



"Loop" secures to bottom of pad. "Hook" to imaging table.

### PERSONALIZATION

#### EMBROIDERY

**#PAD-EML | \$25.00** per pad

Add embroidery to the bottom of your pad. Help identify room numbers, location or even your company logo!\*

[1 line, up to 24 characters available]

**Not available on PAD-51 economy vinyl.**



### PREMIUM PAD STORAGE

#### CORNER STRAPS

**#PAD-ES | \$50.00** per pad

These add-on straps allow hanging of a premium table pad when not in use.

**Not available on PAD-51 economy vinyl or standard pads.**



## Table Pad Storage and Cleaning

Extend the life of your table pad with these accessory options. Storage racks mount on the wall to hang your table pads when not in use. Hanging table pads reduces artifacts and increases the life of your pad!

### TABLE PAD STORAGE OPTIONS

#### STORAGE RACK

**#PAD-R2 | \$63.00** [rack only]

Storage racks easily mount on the wall to hang your table pads when not in use. Hanging your pad reduces wrinkles and increases the life span of your table pad. Ships at 4 lbs. White finish.



Two hooks on 12" centers.  
Mounts on 16" centers.

#### STORAGE RACK & PAD COMBO

**#PAD-S3 | \$265.00** [1 rack & pad combo]

This bundle includes one of our table pad storage racks and a Standard PAD-55 Comfort Foam table pad with black upholstery vinyl. Table pad measures 72" x 23.25" x 1".



Two hooks on 12" centers.  
Mounts on 16" centers.

Easily mounts on the wall to properly store table pads!



### CLEANING AND SANITIZATION

#### TABLE PAD SAFE SANITIZING WIPES

Techno-Aide is proud to offer premium **Sani-Cloth® AF3 Alcohol-Free Wipes**. These disinfectant wipes are perfect for any table pad! These wipes clean, disinfect, deodorize & are fragrance-free. Perfect for other daily use items such as: Sandbags, coated sponges, transfer boards, aprons, shields, vinyl straps and much more!

**Alcohol content:** ALCOHOL-FREE (less than 0.016%). Recommended 3-minutes or less contact Time. Meets CDC, OSHA and CMS Tag F441 guidelines.

CAT. #	Wipe Size	Quantity	Price
GWP-WAF	6" x 6.75"	Tub of 160 wipes	\$20.00
GWP-WAFCS	6" x 6.75"	Sase of 12 tubs	\$235.00



## Table Pad Fabric Options Choose from our unique fabric categories

### UPHOLSTERY VINYL

Durable, Heavy Gauge, Medical-Grade Vinyl, Reinforced 1 mil.



Black  
Upholstery



Green  
Upholstery



Red  
Upholstery



Blue  
Upholstery



Yellow  
Upholstery



Light Blue  
(econ vinyl)

### ECONO-VINYL

Lightweight, Economical Vinyl

### NYLON FABRICS

#### FEATHER RIPSTOP FABRIC Light Weight Nylon, Filament Reinforced Fabric



R701-Royal Blue



R702-Red



R703-Forest



R704-Black



N921-Ocean



N911-Stratos



R706-Navy



R707-Burgundy



R708-Kelly



R709-Charcoal



R710-Gray



N933-Midnight

#### DIAMOND RIPSTOP FABRIC Medium Weight Nylon, Filament Reinforced Fabric



61-Sapphire



64-Onyx



65-Amethyst



62-Ruby



63-Emerald

**NEW ANTI-MICROBIAL FABRIC!**

### TECHNO-TUFF VINYL

Heavyweight, Mesh Laminate, Reinforced Vinyl



TT Black




TT Blue

**Extreme  
Durability &  
Anti-Microbial!**

### HEALTHCARE 4-WAY STRETCH

Comfortable, Pliable, Resilient, Antibacterial & Antifungal



Black-4 way stretch

**Available on premium pads  
and with boxed-ends!**



## Custom Table Pad Form

7117 Centennial Blvd. Nashville, TN 37209

**P (800)251-2629 | P (615)350-7030**

Send all quotes to **Orders@Techno-Aide.com**

TA Part Number

Contact Name	
Address	
City/State/ZIP	
Phone Number	
Email	



Please note, custom orders may take up to 30 business days to manufacture.

A blank sheet of graph paper featuring a uniform grid of small squares. The grid consists of 20 columns and 10 rows, providing a structured space for drawing or writing.

**Pad Style:**

**Foam Type:** ☐ Comfort Foam **Desired Finish:** ☐ Tapered End  
☐ Firm/Bariatric Foam ☐ Tapered w/ Grommets  
☐ Premium Memory Foam ☐ Premium Boxed Ends

<b>Table Pad</b>	<input type="checkbox"/> Economy Vinyl (light blue)	<input type="checkbox"/> Nylon Fabric*	<input type="checkbox"/> Healthcare 4-Way Stretch
<b>Cover Type:</b>	<input type="checkbox"/> Upholstery Vinyl*	<input type="checkbox"/> Techno-Tuff Vinyl*	[Only available in black]

**Color:** \_\_\_\_\_ \*These fabrics have color options PG 11 of our Table Pad PDF

Please write your requested dimensions directly on your sketch and label your lines and angles with the letters below. Please enter the associated dimension into the middle column for each letter used.

Dimension A		Inches
Dimension B		Inches
Dimension C		Inches
Dimension D		Inches
Dimension E		Inches
Dimension F		Inches
Dimension G		Inches
Dimension H		Inches

Economy vinyl cannot have boxed ends or any additional add-ons.  
Only available in light blue.

## TABLE PAD ADD-ONS

- ☐ Add Anti-Slip Backing      ☐ Add Hanging Straps

- ☐ **Add Hook & Loop Strips to the underside of the Pad**

- Hook & Loop **cannot** be combined with "Anti-Slip Backing."
- If "YES," write down whether you want the Hook or the Loop sewn to the pad (This should be the opposite of what is on your table).
- Please note on the diagram above the exact placement of the Hook or Loop and the distances from the sides and ends of the pad, and the length of the strips.

Note: **Hook** is the soft & fuzzy part; **Loop** is the rough & prickly part.

**Customer Signature:** \_\_\_\_\_

Date: \_\_\_\_\_

☐ I confirm that all provided information (including all measurements and attributes) is accurate. I approve manufacturing to these specifications.

## CT Table Slickers

Clear poly film micromatte vinyl protective cover that is designed to fit over CT table pads that come with your table. They are manufactured to resist contamination of the cradle and interior of the imaging table and allows for a quick clean up of blood and bodily fluids and reduces contaminant buildup in hard to clean areas. Flaps prevent contaminants from getting into other hard to clean areas of a table. Such as: table electronics, drives and hinges.

**These slickers are designed to fit the original pad precisely.** Latex Free and Easy to install. Comfortable for patients. Will not interfere with normal CT table operation. Slickers utilize a flap system that is integrated with table safety belts to prevent contamination of CT tables.

If you are covering a non-original table pad, please provide the pad dimensions before placing your order to make sure the slicker cover will fit your pad properly.


## CT Table Cushions

Designed to enhance infection control practices in the CT room. Reduces cleanup times and equipment maintenance. Foam cushion is permanently welded with thermo-sealed seams inside the protective slicker cover, offering contamination resistance of the cradle and interior of the imaging table. The unique flaps prevent contaminants from getting into other hard to clean areas of a table such as threaded surfaces, corners, table joints or cradles, internal electronics, drives and hinges. Cushions are fully enclosed. No need to worry about fitting an existing table pad. Designed to fit the original equipment table specified by machine, these slicker cushions eliminate the fit problem. The most inexpensive way to renew the table pad and cover at once. Allows for a quick clean up of blood and bodily fluids. Reduces contaminant buildup in hard to clean areas. Fluids and contrast are unable to penetrate into the cushion and table. Made of heavy clear micromatte vinyl and top grade foam, fit perfectly to your specific CT table. Easy to install and does not affect normal table usage.

## Features

- Improves table sanitation allowing quicker cleanup
- Reduces contaminant build-up in hard to clean areas
- Reduces table maintenance over time
- Increases system uptime and productivity
- Built to fit original equipment pads precisely!
- Available for nearly every make and model!

CT Table Cushions  
with CT Table Slickers



# Table Pads

CT SLICKERS AND CUSHIONS

TA CAT#	PRICE	SLICKER COVER, SLICKER CUSHION, & ACCESSORY ITEM DESCRIPTION	COMPONENTS	Dimensions (inches)	
GE Slicker Table Covers				L	W
PAD-CT-GE1	\$605	X-ray Angio Slicker Cover for GE Omega IV Tables; fits 3" & 4" thick pad	COV	118"	8"
PAD-CT-GE2	\$605	X-ray Angio/Cardiac Slicker Cover for GE Omega V Tables; fits 3" & 4" thick pad	COV	132"	8"
PAD-CT-GE3	\$605	Angio Slicker - GE Omega IV Tables (118 in.); fits 1" & 2" thick pad	COV V	118"	8"
PAD-CT-GE4	\$605	Angio/Cardiac Slicker - GE Omega V Tables (132 in.); fits 1" & 2" thick pad fits pad E8015JB 1" thick	COV V	132"	8"
PAD-CT-GE5	\$605	CT Slicker w/ Extender - 9800/HiLight Advantage Systems w/ Metal-Free Head Holder	COV E H	78"	16"
PAD-CT-GE6	\$605	CT Slicker w/ Extender - 9800/HiLight Advantage Systems w/o Metal-Free Head Holder	COV E H	75"	16"
PAD-CT-GE7	\$605	CT Slicker - CT/e and CT/e Dual Systems	COV H	82 1/2"	12 1/2"
PAD-CT-GE8	\$605	CT Slicker w/ Extender - BrightSpeed Select & HiSpeed DXi/FXi/LXi/NXi/ZXi Systems	COV E H	86 3/4"	16"
PAD-CT-GE9	\$605	CT Slicker w/ Extender - CT/i & HiSpeed Advantage Systems w/ Metal-Free Head Holder (concave head end)	COV E H	92"	17"
PAD-CT-GE10	\$605	CT Slicker w/ Extender - HiSpeed Advantage Systems w/ Original Head Holder (Square head end)	COV E H	88 1/2"	16"
PAD-CT-GE12	\$605	CT Slicker - ProSpeed Plus Systems	COV G H	79 1/2"	13"
GE Slicker Table Cushions				L	W
PAD-CT-GE15	\$769	CT Slicker Cushion w/ Extender - Discovery CT750 HD and VCT Systems w/ GT 2000 Table	CSH E H	103"	17"
PAD-CT-GE16	\$754	CT Slicker Cushion w/ Extender - Discovery CT750 HD and VCT Systems w/ GT 1700 Table (Optima 660)	CSH E H	82 3/4"	17"
PAD-CT-GE17	\$208	Extender Slicker Pad for VCT, QXi, and PET	E	15"	17"
PAD-CT-GE18	\$711	CT Slicker Cushion w/ Extender - LightSpeed and BrightSpeed Power Systems	CSH E H	92"	17"
PAD-CT-GE20	\$857	PET/CT Slicker Cushion w/ Extender - PET Discovery VCT and Discovery 600/690 (PET Global Table)	CSH E H	110"	17"
PAD-CT-GE21	\$857	PET/CT Slicker Cushion w/ Extender - PET Discovery ST/LS	CSH E H	105"	17"
PAD-CT-GE24	\$711	GE Infinia Nuclear Med Table Cushion	CSH R	77 1/2"	14"
PAD-CT-GE25	\$711	GE Hawkeye Millinium VG Nuclear Med Table Cushion	CSH R	79 5/8"	14"
PAD-CT-GE26	\$711	GE Ventri Nuclear Med Table Cushion (concave head end) fits GE MPR	CSH E R	64"	16"
PAD-CT-GE34	\$766	GE Slicker Cushion w/ Extender - Revolution CT Systems	CSH E H	85 1/2"	17"
PAD-CT-GE41	\$195	GE slicker cushion head pad for CT Optima head section	CSH E H	15"	17"
GE Accessories				L	W
PAD-CT-GE29	\$123	CT Footswitch Slicker - Discovery CT750 HD & VCT Systems (bell shaped)	FSW V	23 1/2"	28"
PAD-CT-GE27	\$73	CT Footswitch Slicker - HiSpeed and LightSpeed Advantage Systems (square shaped)	FSW V	23"	21"
PAD-CT-GE35	\$82	CT Footswitch Slicker - Revolution CT (NG Table) (rectangle shaped)	FSW V	35"	22"
PAD-CT-GE30	\$82	CT Skirt - HiSpeed Advantage Systems (Set of 2 ea.)	S V	13"	36"
Philips Slicker Table Covers				L	W
PAD-CT-PH12	\$605	Philips Nuclear NMPT Pallet Pad Slicker	COV	82 1/2"	15"
Philips Slicker Table Cushions				L	W
PAD-CT-PH13	\$711	Philips Stargate? Mx8000 Dual v.EXP CT Slicker Cushion (Gray)	CSH	94 1/2"	18"
PAD-CT-PH1	\$711	Philips Brilliance CT 40 Channel, 16 Power, 16/10/6 Slice Slicker Cushion (Yellow) Scan line @	CSH V	94 1/2"	18"
PAD-CT-PH2	\$711	Philips Stargate L CT Slicker Cushion Scan line @ 43.31" (111 cm.) from head end	CSH V	94 1/2"	18"
PAD-CT-PH3	\$772	Philips CT Slicker Cushion Long (Yellow) Scan line @ 57.6" (146.3 cm.) from head end	CSH V	104"	18"
PAD-CT-PH4	\$772	Philips CT Slicker Cushion Long (Yellow) Scan line @ 55.5" (141.0 cm.) (Inguinity scanner)	CSH V	104"	18"
PAD-CT-PH5	\$769	Philips Brilliance CT 16/10 Slice CT Slicker Cushion (Yellow) w/ Restraining Belts	CSH V R	94 1/2"	18"
PAD-CT-PH7	\$711	Philips Bariatric Slicker Cushion (Yellow) no flaps	CSH	95"	21"
PAD-CT-PH8	\$772	Philips Bariatric Slicker Cushion, Long (Yellow) no flaps	CSH V	105"	19"
PAD-CT-PH9	\$711	Philips Patient Cushion Scan line at 48.58 inches from head end (Yellow)	CSH V	94 1/2"	18"
PAD-CT-PH10	\$711	Philips Patient Cushion Scan line at 47.24 inches fro head end (Yellow)	CSH V	94 1/2"	18"

PAD-CT-PH14	\$711	Philips Patient Cushion for Flat Top	CSH V	95 1/2"	19
PAD-CT-PH15	\$711	Philips Bariatric Slicker Cushion (Grey)	CSH V	105	19
PAD-CT-PH16	\$813	Mattress Pad-Patient Support Extended	CSH V	106	18"
PAD-CT-PH17	\$789	Mattress Pad-Patient Support Bariatric	CSH V	96 1/2"	18"

**Philips Restraining Belts**

				<b>L</b>	<b>W</b>
PAD-CT-PH11	\$140	Large Restraining belts for Philips Slicker Cushions (Set of 3 ea.)	R	36 1/2"	5 5/8"

**Neusoft Slicker Table Covers**

PAD-CT-NU1	\$605	Slicker cover for Philips and Neusoft	COV		
------------	-------	---------------------------------------	-----	--	--

**Neusoft Slicker Cushions**

PAD-CT-NU2	\$711	Philips and Neusoft CT Slicker Cushion - Grey	CSH		
PAD-CT-NU3	\$711	Philips and Neusoft CT Slicker Cushion - Grey	CSH		
PAD-CT-NU4	\$711	Philips and Neusoft CT Slicker Cushion - Yellow	CSH		
PAD-CT-NU5	\$711	Philips and Neusoft CT Slicker Cushion - Yellow	CSH		

**Siemens Table Slicker Covers**

				<b>L</b>	<b>W</b>
PAD-CT-SM9	\$605	Siemens Somatom Plus 4/Volume Zoom/Sensation CT Slicker Set & Definition DS/AS (HST) MPT Table	COV S H	94 1/2"	18"
PAD-CT-SM6	\$605	Siemens Emotion CT Slicker Set	COV H	86"	18"
PAD-CT-SM7	\$605	Siemens Definition CT Slicker Set	COV H	85"	18"
PAD-CT-SM1	\$605	Siemens Axiom Mattress Shape 1, Model 83 85 333, 4 cm. thick, Angio Slicker	COV	107"	8 1/2"
PAD-CT-SM2	\$605	Siemens Axiom Mattress Shape 1, Model 57 71 444, 8 cm. thick, Angio Slicker	COV	107"	8 1/2"
PAD-CT-SM3	\$605	Siemens Axiom Mattress Shape 2, Model 100 47 902, 4 cm. thick, Angio Slicker	COV	106"	21"
PAD-CT-SM4	\$605	Siemens Axiom Mattress Shape 2, Model 57 71 428, 8 cm. thick, Angio Slicker	COV	106"	21"
PAD-CT-SM30	\$605	Siemens Definition Heavy Duty CT Slicker w/Extender	COV E H	85"	18"
PAD-CT-SM31	\$605	Siemens Sensation Heavy Duty CT Slicker w/Extender	COV S E H	94 1/2"	18"
PAD-CT-SM32	\$605	Siemens Flash CT Slicker Set (Extra wide flaps)	COV H	85"	18"

**Siemens CT Slicker Cushions**

				<b>L</b>	<b>W</b>
PAD-CT-SM16	\$734	Siemens Somatom Plus 4/ VZ/ Sensation/ Definition DS/AS (HST) MPT CT Slicker Cushion	CSH S H	94 1/2"	18"
PAD-CT-SM17	\$711	Siemens Nuclear Medical Slicker Cushion	CSH	83"	14"
PAD-CT-SM18	\$711	Siemens Emotion Slicker Cushion set	CSH H	85 1/2"	18"
PAD-CT-SM20	\$734	Siemens Definition CT Slicker Cushion	CSH H	85"	18"
PAD-CT-SM21	\$1211	Siemens Definition Heavy Duty CT Slicker Cushion Set w/Extender (formally CT-SDEF-06A)	COV PAD E H	85"	18"
PAD-CT-SM22	\$1211	Siemens Sensation Heavy Duty CT Slicker Cushion Set w/Extender	COV PAD S E	94 1/2"	18"

**Siemens Specific Pads**

				<b>L</b>	<b>W</b>
PAD-CT-SM23	\$751	Siemens Definition Heavy Duty CT Flat Top Pad w/Extender (PAD0026XR)	PAD E	85 1/2"	18"
PAD-CT-SM24	\$751	Siemens Sensation Heavy Duty CT Flat Top Pad w/Extender (PAD0026XR)	PAD E	94"	14"

**Siemens Accessories**

				<b>L</b>	<b>W</b>
PAD-CT-SM25	\$208	Siemens Sensation Footboard Extension pad	E	14"	17"
PAD-CT-SM26	\$96	Sterile Drape for Siemens Intervention Panel, 10 drapes per box.	D		
PAD-CT-SM27	\$105	Siemens Somatom Plus 4/VZ/Sensation 54" Skirt Set (Set of 2 ea.)	S	54"	13"
PAD-CT-SM28	\$193	Siemens 64" Restraining Belt (Set of 3 ea.)	R	64"	5 3/4"

## Toshiba CT Slicker Table Covers

				L	W
PAD-CT-TS3	\$605	Toshiba Aquilion CT Slicker (for 81" pad)	COV E H	81"	18"
PAD-CT-TS4	\$605	Toshiba Aquilion CT Slicker (for 70" pad)	COV E H	70"	18"
PAD-CT-TS5	\$605	Toshiba Aquilion CT Slicker Set (for 87" pad)	COV E H	87"	18"
PAD-CT-TS6	\$605	Toshiba Aquilion CT Slicker Set (for 95" pad)	COV E H	95"	18"
PAD-CT-TS7	\$605	Toshiba Astreon CT Slicker (for 75" pad)	COV E H	75"	18"
PAD-CT-TS8	\$605	Toshiba Aquilion One CT Slicker (for 91 1/2" pad)	COV E H	91 1/2"	18"
PAD-CT-TS10	\$605	Slicker Cover for Toshiba Aquilion Prime 96 1/2" Table Pad	COV E H	96 1/2"	19"

## Toshiba Accessories

				L	W
PAD-CT-TS9	\$99	Toshiba Extender Cover	E	15"	18"
PAD-CT-TS11	\$99	Extender cover assembly for Toshiba Aquilion Prime	E	15"	18"

## Miscellaneous

				L	W
PAD-CT-CBH	\$30	All OEMs, 41" Black Catheter Bag Hanger Assembly	H	41"	2"

## Slicker, Cushion and Table Pad Easy Slide Spray

## EASY SLIDE SPRAY - SINGLE

## #CT-ES | \$20 Single Bottle

This water-based silicone spray increases the slickness of a desired surface, reducing the effort required for patient transfer and positioning! Ideal for table pads, slickers and cushions! This spray is safe for skin contact, non-flammable and produces no smell or fumes. Designed to work on all plastic and vinyl surfaces, Easy Slide spray comes in a convenient pump spray bottle. Available as a single bottle or in a case of 12 bottles.



## EASY SLIDE SPRAY - CASE

## #CT-ES12 | \$225 Case of 12 bottles

Easy Slide spray is now available in a case of 12 individual bottles.





ADMINISTRATION/MANUFACTURING

7117 Centennial Blvd. Nashville, TN 37209

Orders can be placed online or sent to [orders@techno-aide.com](mailto:orders@techno-aide.com)

For official Quotes: [quotes@techno-aide.com](mailto:quotes@techno-aide.com)

Any questions can be directed to [customerservice@techno-aide.com](mailto:customerservice@techno-aide.com)

P: 800.251.2629 or 615.350.7030

F: 800.554.6275 or 615.350.7879

[techno-aide.com](http://techno-aide.com)