

X-ray Generator



System Model	VXR-32	VXR-40	VXR-C32	VXR-C40
Generator Model	GXR-32	GXR-40	GXR-C32	GXR-C40
Output Rating	32kW	40kW	32kW	40kW
Line Power Rating	230VAC, Single phase, ±10% (50*/60Hz), * Outside North America		110~120, 220~230VAC, Single phase, ±10% (50*/60Hz), 3kVA * Outside North America	
kVp Range	40~125kVp, 1kV step		40~125kVp, 1kV step	
mA Range	10 to 400mA	10 to 500mA	10 to 400mA	10 to 500mA
mAs Range	0.1 ~ 500mAs		0.1 ~ 500mAs	
Time Range	1msec ~ 10 sec		1msec ~ 10 sec	

X-ray Tube Assembly


Focal Spot Size	1.0 / 2.0 mm	0.6 / 1.2 mm	0.6/1.5 mm
Max. Tube Voltage	125 kV	125 kV	125 kV
Max. Anode HU	140kHU	150kHU	200 kHU

Tube Stand & Table

Tabletop Size	1,400(W) x 700(D) x 47(H)mm
Tabletop Longitudinal Travel	400(±200)mm
SID	1,100mm
Dimension	1,400(W) x 984(D) x 1,940(H)mm

Dimension

Unit : mm (inch)



Simple & Easy Operation !
Excellent Image Quality !!



VXR SERIES
Veterinary X-ray System