


TechnoAide
THE ONE-STOP-SHOP FOR IMAGING SOLUTIONS



Catalog '24

Catalog '24

Catalog '24

Catalog '24

Catalog '24

Catalog '24

Table of Contents



Protective Apparel & Supplies **1**



Sponges & Positioning Aids **27**



Veterinary Positioning Aids **40**



Table Pads **41**



X-Ray Markers **43**






General Supplies & Accessories **45**



Index **47**

MEASURING MADE EZ.



Measuring lead shouldn't be complicated, it should be EZ!
Try our new EZ Fit Apparel Sizing Tool

FREE!

Scan or visit TECHNO-AIDE.com to learn more.



APRON·HUB

Techno-Aide's revolutionary **FREE** lead tracking system. Now with built-in SmartID® data import capability!

WEB-BASED ACCESS

from any device.

UNIVERSALLY COMPATIBLE

with any QR/Barcode Scanner.

ACCOUNT CUSTOMIZATION TOOLS

to meet the specific needs of your department.

CUSTOM REPORTING

makes exporting inventory data specific and easy.

MASS IMPORTING OPTION

to make switching from any platform a streamlines process.


DIGITAL RETAGGING

features a first-time tool built right into each product page to allow for seamless apron retagging.

CLOUD-BASED SECURITY

for built in scale no matter how big your organization becomes.

ABSOLUTELY FREE!



ApronHub.com
CustomerService@Techno-Aide.com
1-800-251-2629



EZ FRONT APRON

Lead Options:

- Regular
- Lightweight
- Super-Lite® Lead-Free

Available in:

- 0.35mm Pb front protection
- 0.5mm Pb front protection

Features:

- Crisscross strap design & shoulder strap for increased comfort
- Unisex Sizing
- Precision-cut necklines & premium shoulder pads for comfort & durability
- Hook & Loop or Buckle Closure (no additional cost)
- Optional embroidery/monogramming (no additional cost; see pg. 23 for details)
- Optional Thyroid buckle tether



| 0.35mm | Regular | Lightweight | Super-Lite® Lead-Free |
|---------|---------|-------------|-----------------------|
| Petite | EP-3R | EP-3L | EP-3N |
| Small | ES-3R | ES-3L | ES-3N |
| Medium | EM-3R | EM-3L | EM-3N |
| Large | EL-3R | EL-3L | EL-3N |
| X-Large | EX-3R | EX-3L | EX-3N |
| 2XL | E2-3R | E2-3L | E2-3N |
| 3XL | E3-3R | E3-3L | E3-3N |
| 0.5mm | Regular | Lightweight | Super-Lite® Lead-Free |
| Petite | EP-5R | EP-5L | EP-5N |
| Small | ES-5R | ES-5L | ES-5N |
| Medium | EM-5R | EM-5L | EM-5N |
| Large | EL-5R | EL-5L | EL-5N |
| X-Large | EX-5R | EX-5L | EX-5N |
| 2XL | E2-5R | E2-5L | E2-5N |
| 3XL | E3-5R | E3-5L | E3-5N |



EZ COMFORT-FLEX FRONT APRON

Lead Options:

- Regular
- Lightweight
- Super-Lite® Lead-Free

Available in:

- 0.35mm Pb front protection
- 0.5mm Pb front protection

Features:

- Superior flexible crisscross strap design & shoulder strap for increased comfort
- Unisex Sizing
- Precision-cut necklines & premium shoulder pads for comfort & durability
- Hook & Loop or Buckle Closure (no additional cost)
- Optional embroidery/monogramming (no additional cost; see pg. 23 for details)
- Optional Thyroid buckle tether



| 0.35mm | Regular | Lightweight | Super-Lite® Lead-Free |
|---------|---------|-------------|-----------------------|
| Petite | FP-3R | FP-3L | FP-3N |
| Small | FS-3R | FS-3L | FS-3N |
| Medium | FM-3R | FM-3L | FM-3N |
| Large | FL-3R | FL-3L | FL-3N |
| X-Large | FX-3R | FX-3L | FX-3N |
| 2XL | F2-3R | F2-3L | F2-3N |
| 3XL | F3-3R | F3-3L | F3-3N |
| 0.5mm | Regular | Lightweight | Super-Lite® Lead-Free |
| Petite | FP-5R | FP-5L | FP-5N |
| Small | FS-5R | FS-5L | FS-5N |
| Medium | FM-5R | FM-5L | FM-5N |
| Large | FL-5R | FL-5L | FL-5N |
| X-Large | FX-5R | FX-5L | FX-5N |
| 2XL | F2-5R | F2-5L | F2-5N |
| 3XL | F3-5R | F3-5L | F3-5N |



EZ RELEASE SURGERY FRONT APRON

Lead Options:

- Regular
- Lightweight
- Super-Lite® Lead-Free

Available in:

- 0.35mm Pb front protection
- 0.5mm Pb front protection

Features:

- Quick-drop cross-strap design for optional hands-free removal
- Unisex Sizing
- Precision-cut necklines for comfort & durability
- Optional embroidery/monogramming (no additional cost; see pg. 23 for details)
- Optional Thyroid buckle tether



| 0.35mm | Regular | Lightweight | Super-Lite® Lead-Free |
|---------|---------|-------------|-----------------------|
| Petite | UP-3R | UP-3L | UP-3N |
| Small | US-3R | US-3L | US-3N |
| Medium | UM-3R | UM-3L | UM-3N |
| Large | UL-3R | UL-3L | UL-3N |
| X-Large | UX-3R | UX-3L | UX-3N |
| 2XL | U2-3R | U2-3L | U2-3N |
| 3XL | U3-3R | U3-3L | U3-3N |
| 0.5mm | Regular | Lightweight | Super-Lite® Lead-Free |
| Petite | UP-5R | UP-5L | UP-5N |
| Small | US-5R | US-5L | US-5N |
| Medium | UM-5R | UM-5L | UM-5N |
| Large | UL-5R | UL-5L | UL-5N |
| X-Large | UX-5R | UX-5L | UX-5N |
| 2XL | U2-5R | U2-5L | U2-5N |
| 3XL | U3-5R | U3-5L | U3-5N |



EZ RELEASE COMFORT OR FRONT APRON

Lead Options:

- Regular
- Lightweight
- Super-Lite® Lead-Free

Available in:

- 0.35mm Pb front protection
- 0.5mm Pb front protection

Features:

- Dual quick-release shoulder straps for shoulder relief and optional hands-free removal
- Unisex Sizing
- Precision-cut necklines & premium shoulder pads for comfort & durability
- Optional embroidery/monogramming (no additional cost; see pg. 23 for details)
- Optional Thyroid buckle tether



| 0.35mm | Regular | Lightweight | Super-Lite® Lead-Free |
|---------|---------|-------------|-----------------------|
| Petite | OP-3R | OP-3L | OP-3N |
| Small | OS-3R | OS-3L | OS-3N |
| Medium | OM-3R | OM-3L | OM-3N |
| Large | OL-3R | OL-3L | OL-3N |
| X-Large | OX-3R | OX-3L | OX-3N |
| 2XL | O2-3R | O2-3L | O2-3N |
| 3XL | O3-3R | O3-3L | O3-3N |
| 0.5mm | Regular | Lightweight | Super-Lite® Lead-Free |
| Petite | OP-5R | OP-5L | OP-5N |
| Small | OS-5R | OS-5L | OS-5N |
| Medium | OM-5R | OM-5L | OM-5N |
| Large | OL-5R | OL-5L | OL-5N |
| X-Large | OX-5R | OX-5L | OX-5N |
| 2XL | O2-5R | O2-5L | O2-5N |
| 3XL | O3-5R | O3-5L | O3-5N |



EZ FULL WRAP VEST

Lead Options:

- Regular
- Lightweight
- Super-Lite® Lead-Free

Available in:

- 0.35mm Pb front protection & 0.175mm Pb back protection
- 0.5mm Pb front protection & 0.25mm Pb back protection

Features:

- Full panel overlap for complete frontal protection
- Unisex Sizing
- Precision-cut necklines & premium shoulder pads for comfort & durability
- Hook & Loop or Combo Closure (no additional cost)
- Optional embroidery/monogramming (no additional cost; see pg. 23 for details)
- Optional Thyroid buckle tether



| 0.35mm | Regular | Lightweight | Super-Lite® Lead-Free |
|---------|---------|-------------|-----------------------|
| Petite | VP-3R | VP-3L | VP-3N |
| Small | VS-3R | VS-3L | VS-3N |
| Medium | VM-3R | VM-3L | VM-3N |
| Large | VL-3R | VL-3L | VL-3N |
| X-Large | VX-3R | VX-3L | VX-3N |
| 2XL | V2-3R | V2-3L | V2-3N |
| 3XL | V3-3R | V3-3L | V3-3N |
| 0.5mm | Regular | Lightweight | Super-Lite® Lead-Free |
| Petite | VP-5R | VP-5L | VP-5N |
| Small | VS-5R | VS-5L | VS-5N |
| Medium | VM-5R | VM-5L | VM-5N |
| Large | VL-5R | VL-5L | VL-5N |
| X-Large | VX-5R | VX-5L | VX-5N |
| 2XL | V2-5R | V2-5L | V2-5N |
| 3XL | V3-5R | V3-5L | V3-5N |



EZ FULL WRAP CONTOUR VEST

Lead Options:

- Regular
- Lightweight
- Super-Lite® Lead-Free

Available in:

- 0.35mm Pb front protection & 0.175mm Pb back protection
- 0.5mm Pb front protection & 0.25mm Pb back protection

Features:

- Improved chest coverage, specifically for women
- Unisex Sizing
- Tapered cut for improved weight distribution & comfort
- Less gapping around shoulders & neck
- Hook & Loop or Combo Closure (no additional cost)
- Optional embroidery/monogramming (no additional cost; see pg. 23 for details)
- Optional Thyroid buckle tether



| 0.35mm | Regular | Lightweight | Super-Lite® Lead-Free |
|---------|---------|-------------|-----------------------|
| Petite | XP-3R | XP-3L | XP-3N |
| Small | XS-3R | XS-3L | XS-3N |
| Medium | XM-3R | XM-3L | XM-3N |
| Large | XL-3R | XL-3L | XL-3N |
| X-Large | XX-3R | XX-3L | XX-3N |
| 2XL | X2-3R | X2-3L | X2-3N |
| 3XL | X3-3R | X3-3L | X3-3N |
| 0.5mm | Regular | Lightweight | Super-Lite® Lead-Free |
| Petite | XP-5R | XP-5L | XP-5N |
| Small | XS-5R | XS-5L | XS-5N |
| Medium | XM-5R | XM-5L | XM-5N |
| Large | XL-5R | XL-5L | XL-5N |
| X-Large | XX-5R | XX-5L | XX-5N |
| 2XL | X2-5R | X2-5L | X2-5N |
| 3XL | X3-5R | X3-5L | X3-5N |



EZ REVERSE FULL WRAP VEST

Lead Options:

- Regular
- Lightweight
- Super-Lite® Lead-Free

Available in:

- 0.35mm Pb front & back protection
- 0.5mm Pb front & back protection

Features:

- Unique reverse wrap design for complete front & back protection
- Unisex Sizing
- Precision-cut necklines & premium shoulder pads for comfort & durability
- Optional embroidery/monogramming (no additional cost; see pg. 23 for details)
- Optional Thyroid buckle tether



| 0.35mm | Regular | Lightweight | Super-Lite® Lead-Free |
|---------|---------|-------------|-----------------------|
| Petite | WP-3R | WP-3L | WP-3N |
| Small | WS-3R | WS-3L | WS-3N |
| Medium | WM-3R | WM-3L | WM-3N |
| Large | WL-3R | WL-3L | WL-3N |
| X-Large | WX-3R | WX-3L | WX-3N |
| 2XL | W2-3R | W2-3L | W2-3N |
| 3XL | W3-3R | W3-3L | W3-3N |
| 0.5mm | Regular | Lightweight | Super-Lite® Lead-Free |
| Petite | WP-5R | WP-5L | WP-5N |
| Small | WS-5R | WS-5L | WS-5N |
| Medium | WM-5R | WM-5L | WM-5N |
| Large | WL-5R | WL-5L | WL-5N |
| X-Large | WX-5R | WX-5L | WX-5N |
| 2XL | W2-5R | W2-5L | W2-5N |
| 3XL | W3-5R | W3-5L | W3-5N |



EZ FULL WRAP SKIRT

Lead Options:

- Regular
- Lightweight
- Super-Lite® Lead-Free

Available in:

- 0.35mm Pb front protection & 0.175mm Pb back protection
- 0.5mm Pb front protection & 0.25mm Pb back protection

Features:

- Full panel overlap for complete frontal protection
- Even weight distribution across waist and hips for extended comfort
- Unisex Sizing
- Combo Closure (no additional cost)
- Optional embroidery/monogramming (no additional cost; see pg. 23 for details)



| 0.35mm | Regular | Lightweight | Super-Lite® Lead-Free |
|---------|---------|-------------|-----------------------|
| Petite | SP-3R | SP-3L | SP-3N |
| Small | SS-3R | SS-3L | SS-3N |
| Medium | SM-3R | SM-3L | SM-3N |
| Large | SL-3R | SL-3L | SL-3N |
| X-Large | SX-3R | SX-3L | SX-3N |
| 2XL | S2-3R | S2-3L | S2-3N |
| 3XL | S3-3R | S3-3L | S3-3N |
| 0.5mm | Regular | Lightweight | Super-Lite® Lead-Free |
| Petite | SP-5R | SP-5L | SP-5N |
| Small | SS-5R | SS-5L | SS-5N |
| Medium | SM-5R | SM-5L | SM-5N |
| Large | SL-5R | SL-5L | SL-5N |
| X-Large | SX-5R | SX-5L | SX-5N |
| 2XL | S2-5R | S2-5L | S2-5N |
| 3XL | S3-5R | S3-5L | S3-5N |



EZ FULL WRAP APRON

Lead Options:

- Lightweight
- Super-Lite® Lead-Free

Available in:

- 0.35mm Pb front protection & 0.175mm Pb back protection
- 0.5mm Pb front protection & 0.25mm Pb back protection

Features:

- Full panel overlap for complete frontal protection
- Combo closure (hook & loop and buckle) design for even weight distribution across shoulders, waist, and hips
- Unisex Sizing
- Precision-cut necklines & premium shoulder pads for comfort & durability
- Combo Closure (no additional cost)
- Optional embroidery/monogramming (no additional cost; see pg. 23 for details)
- Optional Thyroid buckle tether



| 0.35mm | Lightweight | Super-Lite® Lead-Free |
|---------|-------------|-----------------------|
| Petite | MP-3L | MP-3N |
| Small | MS-3L | MS-3N |
| Medium | MM-3L | MM-3N |
| Large | ML-3L | ML-3N |
| X-Large | MX-3L | MX-3N |
| 2XL | M2-3L | M2-3N |
| 3XL | M3-3L | M3-3N |
| 0.5mm | Lightweight | Super-Lite® Lead-Free |
| Petite | MP-5L | MP-5N |
| Small | MS-5L | MS-5N |
| Medium | MM-5L | MM-5N |
| Large | ML-5L | ML-5N |
| X-Large | MX-5L | MX-5N |
| 2XL | M2-5L | M2-5N |
| 3XL | M3-5L | M3-5N |



Thyroid Protection

Techno-Aide's **Standard Thyroid Guards** are a convenient and cost-effective solution for thyroid protection. Thyroid Guards feature a simple one-size-fits-all design and precision-cut comfort necklines. Available in three styles to meet virtually any need: top and bottom "Full Binding", bottom only "Mid-Binding", or premium "No Binding." All Thyroid Guards are available with either **0.5mm Pb** or **0.35mm Pb equivalent protection**, and your choice of three media options: Regular Lead, Lightweight Lead (composite) or our premium Super-Lite® Lead-Free. All Standard Thyroid Guards include free embroidery/monogramming (see pg. 23 for details) and a convenient buckle tether.

— Choose from 3 **FREE** closure options: Buckle, Hook & Loop, or Magnetic —

FULL BINDING

| 0.35mm | Regular | Lightweight | Super-Lite® Lead-Free |
|-------------|---------|-------------|-----------------------|
| Buckle | TS-3R-B | TS-3L-B | TS-3N-B |
| Hook & Loop | TS-3R-V | TS-3L-V | TS-3N-V |
| Magnetic | TS-3R-M | TS-3L-M | TS-3N-M |
| 0.5mm | Regular | Lightweight | Super-Lite® Lead-Free |
| Buckle | TS-5R-B | TS-5L-B | TS-5N-B |
| Hook & Loop | TS-5R-V | TS-5L-V | TS-5N-V |
| Magnetic | TS-5R-M | TS-5L-M | TS-5N-M |



MID-BINDING

| 0.35mm | Regular | Lightweight | Super-Lite® Lead-Free |
|-------------|---------|-------------|-----------------------|
| Buckle | TH-3R-B | TH-3L-B | TH-3N-B |
| Hook & Loop | TH-3R-V | TH-3L-V | TH-3N-V |
| Magnetic | TH-3R-M | TH-3L-M | TH-3N-M |
| 0.5mm | Regular | Lightweight | Super-Lite® Lead-Free |
| Buckle | TH-5R-B | TH-5L-B | TH-5N-B |
| Hook & Loop | TH-5R-V | TH-5L-V | TH-5N-V |
| Magnetic | TH-5R-M | TH-5L-M | TH-5N-M |



NO BINDING

| 0.35mm | Regular | Lightweight | Super-Lite® Lead-Free |
|-------------|---------|-------------|-----------------------|
| Buckle | TN-3R-B | TN-3L-B | TN-3N-B |
| Hook & Loop | TN-3R-V | TN-3L-V | TN-3N-V |
| Magnetic | TN-3R-M | TN-3L-M | TN-3N-M |
| 0.5mm | Regular | Lightweight | Super-Lite® Lead-Free |
| Buckle | TN-5R-B | TN-5L-B | TN-5N-B |
| Hook & Loop | TN-5R-V | TN-5L-V | TN-5N-V |
| Magnetic | TN-5R-M | TN-5L-M | TN-5N-M |





Techno-Aide is proud to offer two **Specialty Thyroid Guards** styles to help meet unique patient needs. These Specialty Thyroid Guard styles feature a simple one-size-fits-all design and precision-cut comfort necklines. Specialty Thyroid Guards feature **0.5mm Pb equivalent protection**, and your choice of three media options: Regular Lead, Lightweight Lead (composite) or our premium Super-Lite® Lead-Free. All Specialty Thyroid Guards include free embroidery/monogramming (see pg. 23 for details).

MAMMOGRAPHY

Techno-Aide's **Mammography Thyroid Guards** feature a premium "No Binding" style for maximum comfort. Designed for use in mammography studies, the short bib protects the patient's thyroid without interfering with the mammography image.

| 0.5mm | Regular | Lightweight | Super-Lite® Lead-Free |
|-------|---------|-------------|-----------------------|
| Mammo | TM-5R | TM-5L | TM-5N |

PEDIATRIC

Techno-Aide's **Pediatric Thyroid Guards** are designed for use in pediatric studies. The one-size-fits-all design and compact bib are perfect for adolescent patients of virtually any age.

| 0.5mm | Regular | Lightweight | Super-Lite® Lead-Free |
|-----------|---------|-------------|-----------------------|
| Pediatric | TP-5R | TP-5L | TP-5N |



All **FULL, MID, and NO BINDING** Thyroid Guards include the option of a **FREE** buckle tether!



Techno-Aide is proud to offer our unique line of **General and Specialty Protective Apparel**. This diverse line of products and accessories are specially curated to help meet a variety of unique patient and facility needs.

LAP GUARD HALF APRONS

| 0.35mm | Regular | Lightweight | Super-Lite® Lead-Free |
|---------|---------|-------------|--------------------------|
| Petite | LP-3R | LP-3L | LP-3N |
| Small | LS-3R | LS-3L | LS-3N |
| Medium | LM-3R | LM-3L | LM-3N |
| Large | LL-3R | LL-3L | LL-3N |
| X-Large | LX-3R | LX-3L | LX-3N |
| 2XL | L2-3R | L2-3L | L2-3N |
| 0.5mm | Regular | Lightweight | Super-Lite® Lead-Free |
| Petite | LP-5R | LP-5L | LP-5N |
| Small | LS-5R | LS-5L | LS-5N |
| Medium | LM-5R | LM-5L | LM-5N |
| Large | LL-5R | LL-5L | LL-5N |
| X-Large | LX-5R | LX-5L | LX-5N |
| 2XL | L2-5R | L2-5L | L2-5N |




Techno-Aide is proud to offer our **Lap Guard Half Apron Bundles**. These prebuilt bundles are a convenient and quick way to outfit virtually any imaging facility. Each bundle includes four color-coded lap guard half apron sizes and a low-profile apron storage rack. Available with your choice of hook & loop or buckle closures and free embroidery.

LAP GUARD HALF APRON BUNDLES

Included Products:

- S** Lap Guard Half Apron - **Berry** - 0.5mm Pb Reg. Lead
- M** Lap Guard Half Apron - **Gray** - 0.5mm Pb Reg. Lead
- L** Lap Guard Half Apron - **Royal Blue** - 0.5mm Pb Reg. Lead
- XL** Lap Guard Half Apron - **Black** - 0.5mm Pb Reg. Lead
- Wall Mounted Storage Rack** - 2 Pegs [white]

| SKU | Buckle | Hook & Loop | Embroidery |
|--------------|--------|-------------|------------|
| LAP-BNDL-B | ✓ | | |
| LAP-BNDL-B-E | ✓ | | ✓ |
| LAP-BNDL-V | | ✓ | |
| LAP-BNDL-V-E | | ✓ | ✓ |





RADIATION SAFETY GLASSES



NIKE brand available!

– Now offering prescription lenses! –

| SKU | Part Description |
|-------------------|---|
| RG-208-BK | Ultra-Guard Radiation Safety Glasses - Plano Lens (Black) |
| RG-208-S | Ultra-Guard Radiation Safety Glasses - Plano Lens (Silver) |
| RG-208-R | Ultra-Guard Radiation Safety Glasses - Plano Lens (Red) |
| RG-208-TA | Ultra-Guard Radiation Safety Glasses - Plano Lens (Taupe) |
| RG-208-PG | Ultra-Guard Radiation Safety Glasses - Plano Lens (Purple/Gray) |
| RG-208-BY | Ultra-Guard Radiation Safety Glasses - Plano Lens (Blue/Yellow) |
| RG-208-BK | Ultra-Guard Radiation Safety Glasses - Plano Lens (Pink) |
| RG-33-BK-50SS | Cover-Guard Radiation Safety Glasses - Plano Lens (Black) |
| RG-33-T-50SS | Cover-Guard Radiation Safety Glasses - Plano Lens (Tortoise) |
| RG-33-W-50SS | Cover-Guard Radiation Safety Glasses - Plano Lens (White) |
| RG-1205-BK | Rugged Radiation Safety Glasses - Black |
| RG-1205-BL | Rugged Radiation Safety Glasses - Blue |
| RG-15011-BKC-50SS | Premium Radiation Safety Glasses - Black Crystal |
| RG-15011-BKY-50SS | Premium Radiation Safety Glasses - Black Green |
| RG-15011-CR-50SS | Premium Radiation Safety Glasses - Crystal Red |
| RG-15011-CC-50SS | Premium Radiation Safety Glasses - Crystal Clear |
| RG-15011-NVC-50SS | Premium Radiation Safety Glasses - Navy Clear |
| RG-16001-BKBL | Deluxe Radiation Safety Glasses - Black/Blue |
| RG-F10-BB-50SS | Economy Radiation Safety Glasses - Blue Black |



RG-F10-CLR-50SS



RG-1205-BL



RG-NI-FJ1975-010

| SKU | Part Description |
|-------------------|---|
| RG-F10-BY-50SS | Economy Radiation Safety Glasses - Black Yellow |
| RG-F10-CLR-50SS | Economy Radiation Safety Glasses - Black Yellow Clear |
| RG-F126-BK-50SS | Standard Radiation Safety Glasses - Black |
| RG-NI-FJ1994-060 | Nike Brazen Boost Radiation Safety Glasses - Matte Anthracite |
| RG-NI-FJ1978-451 | Nike Brazen Boost Radiation Safety Glasses - Matte Dark Obsidian |
| RG-NI-FJ1975-010 | Nike Brazen Boost Radiation Safety Glasses - Matte Black/White |
| RG-NI-FJ1975-021 | Nike Brazen Boost Radiation Safety Glasses - Matte Dark Gray/White |
| RG-NI-FJ1985-012 | Nike Brazen Shadow Radiation Safety Glasses - Matte Dark Gray/Gray |
| RG-NI-FJ1987-451 | Nike Brazen Shadow Radiation Safety Glasses - Matte Dark Obsidian |
| RG-NI-FJ1985-060 | Nike Brazen Shadow Radiation Safety Glasses - Matte Anthracite Gray |
| RG-NI-FJ1985-010 | Nike Brazen Shadow Radiation Safety Glasses - Matte Black/White |
| RG-NI-FJ1987-021 | Nike Brazen Shadow Radiation Safety Glasses - Matte Dark Gray/Red |
| RG-NI-EV0571-001 | Nike Brazen Radiation Safety Glasses - Black |
| RG-NI-EV0758-049 | Nike Brazen Radiation Safety Glasses - Matte Black/Military Blue |
| RG-Q100-BK-50SS | Convertible Radiation Safety Glasses - Black/Clear [Large] |
| RG-Q200-BK-50SS | Convertible Radiation Safety Glasses - Glossy Black [Medium] |
| RG-Q200-BLPK-50SS | Convertible Radiation Safety Glasses - Blue/Pink [Medium] |
| RG-Q200-PUR-50SS | Convertible Radiation Safety Glasses - Purple [Medium] |
| RG-Q300-BK-50SS | Convertible Radiation Safety Glasses - Matte Black [Small] |

RADIATION SHIELDING GLOVES



100V-DKB



101SLMIT-DKB



305101




305202

| SKU | Part Description |
|--------------|---|
| 100V-DKB | Premium Full Coverage Radiation Safety Finger Gloves (Single Pair) |
| 101SLMIT-DKB | Premium Full Coverage Radiation Safety Mitts w/ Dexterity Palm Slit (Single Pair) |
| 305101 | Lite Surgical Radition Protection Gloves - 0.20mm [Size 6.5 - XS] (QTY 5 Pairs) |
| 305102 | Lite Surgical Radition Protection Gloves - 0.20mm [Size 7 - S] (QTY 5 Pairs) |
| 305103 | Lite Surgical Radition Protection Gloves - 0.20mm [Size 7.5 - M] (QTY 5 Pairs) |
| 305104 | Lite Surgical Radition Protection Gloves - 0.20mm [Size 8 - L] (QTY 5 Pairs) |
| 305105 | Lite Surgical Radition Protection Gloves - 0.20mm [Size 8.5 - XL] (QTY 5 Pairs) |
| 305202 | Max Surgical Radition Protection Gloves - 0.40mm [Size 7 - S] (QTY 5 Pairs) |
| 305203 | Max Surgical Radition Protection Gloves - 0.40mm [Size 7.5 - M] (QTY 5 Pairs) |
| 305204 | Max Surgical Radition Protection Gloves - 0.40mm [Size 8 - L] (QTY 5 Pairs) |



BLANKET GUARDS

| 0.35mm | Regular | Lightweight | Super-Lite® Lead-Free |
|---------|---------|-------------|--------------------------|
| Small | AS-3R | AS-3L | AS-3N |
| Medium | AM-3R | AM-3L | AM-3N |
| Large | AL-3R | AL-3L | AL-3N |
| X-Large | AX-3R | AX-3L | AX-3N |
| 0.5mm | Regular | Lightweight | Super-Lite® Lead-Free |
| Small | AS-5R | AS-5L | AS-5N |
| Medium | AM-5R | AM-5L | AM-5N |
| Large | AL-5R | AL-5L | AL-5N |
| X-Large | AX-5R | AX-5L | AX-5N |



GONADAL GUARDS

| 0.5mm | Regular | Lightweight | Super-Lite® Lead-Free |
|---------|---------------|---------------|--------------------------|
| Small | GS-5R | GS-5L | GS-5N |
| Medium | GM-5R | GM-5L | GM-5N |
| Large | GL-5R | GL-5L | GL-5N |
| Bundles | GB-5R (1 ea.) | GB-5L (1 ea.) | GB-5N (1 ea.) |




BREAST GUARDS

| 0.5mm | Regular | Lightweight | Super-Lite® Lead-Free |
|-------|---------|-------------|--------------------------|
| Child | BC-5R | BC-5L | BC-5N |
| Adult | BA-5R | BA-5L | BA-5N |



SCOLIO GUARDS

| 0.5mm | Regular |
|----------|---------|
| Standard | IS-5R |





PREGNANCY/MATERNITY APRONS

PREGNANCY FRONT APRONS

| 1.0mm | Lightweight | Super-Lite® Lead-Free |
|----------|-------------|--------------------------|
| Standard | PS-1L | PS-1N |
| Full | PT-1L | PT-1N |

PREGNANCY LAP GAUARD HALF APRONS

| 1.0mm | Lightweight | Super-Lite® Lead-Free |
|----------|-------------|--------------------------|
| Standard | RX-1L | RX-1N |



EZ PEDIATRIC FRONT APRON

| 0.5mm | Regular | Lightweight | Super-Lite® Lead-Free |
|-------|---------|-------------|--------------------------|
| Child | CC-3R | CC-3L | CC-3N |



DENTAL APRONS

EZ DENTAL DRAPES

| 0.35mm | Regular | Lightweight | Super-Lite® Lead-Free |
|--------|---------|-------------|--------------------------|
| Child | DC-3R | DC-3L | DC-3N |
| Adult | DA-3R | DA-3L | DA-3N |

EZ PANORAMIC DENTAL DRAPES

| 0.35mm | Regular | Lightweight | Super-Lite® Lead-Free |
|--------|---------|-------------|--------------------------|
| Adult | NA-3R | NA-3L | NA-3N |





MOBILE SHIELDS

| SKU | Part Description |
|------------|--|
| LPS-S2424 | Economy Star-Base Mobile Shield [24 x 24/0.5mm Regular Lead/5-Star Base] |
| LPS-S2436 | Economy Star-Base Mobile Shield [24 x 36/0.5mm Regular Lead/5-Star Base] |
| LPS-S2448 | Economy Star-Base Mobile Shield [24 x 48/0.5mm Regular Lead/5-Star Base] |
| LPSG-S2424 | Economy Star-Base Mobile Shield [24 x 24/0.5mm Super-Lite® Lead-Free/5-Star Base] |
| LPSG-S2436 | Economy Star-Base Mobile Shield [24 x 36/0.5mm Super-Lite® Lead-Free/5-Star Base] |
| LPSG-S2448 | Economy Star-Base Mobile Shield [24 x 48/0.5mm Super-Lite® Lead-Free/5-Star Base] |
| LPS-T2424S | Standard T-Base Mobile Shield [24 x 24/0.5mm Regular Lead/T-Base] |
| LPS-T2430S | Standard T-Base Mobile Shield [24 x 30/0.5mm Regular Lead/T-Base] |
| LPS-T2436S | Standard T-Base Mobile Shield [24 x 36/0.5mm Regular Lead/T-Base] |
| LPS-T2448S | Standard T-Base Mobile Shield [24 x 48/0.5mm Regular Lead/T-Base] |
| LPS-T3030S | Standard T-Base Mobile Shield [30 x 30/0.5mm Regular Lead/T-Base] |
| LPS-T3036S | Standard T-Base Mobile Shield [30 x 36/0.5mm Regular Lead/T-Base] |
| LPS-T3042S | Standard T-Base Mobile Shield [30 x 42/0.5mm Regular Lead/T-Base] |
| LPS-H3636 | Premium H-Base Mobile Shield [36 x 36/0.5mm Regular Lead/H-Base] |
| LPS-H3648 | Premium H-Base Mobile Shield [36 x 48/0.5mm Regular Lead/H-Base] |
| LPS-H4848 | Premium H-Base Mobile Shield [48 x 48/0.5mm Regular Lead/H-Base] |



LEADED TABLE DRAPES

| SKU | Part Description |
|-----------|---|
| TD1-18-11 | Leaded Table Drape - Single 18" x 12" Panel (L x W) 0.5mm Regular Lead 11 - Royal Blue |
| TD1-18-33 | Leaded Table Drape - Single 18" x 12" Panel (L x W) 0.5mm Regular Lead 33 - Black |
| TD1-24-11 | Leaded Table Drape - Single 24" x 12" Panel (L x W) 0.5mm Regular Lead 11 - Royal Blue |
| TD1-24-33 | Leaded Table Drape - Single 24" x 12" Panel (L x W) 0.5mm Regular Lead 33 - Black |
| TD1-30-11 | Leaded Table Drape - Single 30" x 12" Panel (L x W) 0.5mm Regular Lead 11 - Royal Blue |
| TD1-30-33 | Leaded Table Drape - Single 30" x 12" Panel (L x W) 0.5mm Regular Lead 33 - Black |
| TD1-36-11 | Leaded Table Drape - Single 36" x 12" Panel (L x W) 0.5mm Regular Lead 11 - Royal Blue |
| TD1-36-33 | Leaded Table Drape - Single 36" x 12" Panel (L x W) 0.5mm Regular Lead 33 - Black |
| TL-1-24 | Table Drape Adhesive Table Liner - Single 24" x 11" Strip (L x W) [Black] |
| TL-1-48 | Table Drape Adhesive Table Liner - Single 48" x 12" Panel (L x W) [Black] |
| TDHS | Table Drape Hanging Straps - Set of 2 (8" Hook & Loop Straps) [Black] |





STANDARD MOBILE BARRIERS


LB-2430

LB-3048

LB-3060

LB-3036

| SKU | Part Description |
|---------|---|
| LB-2430 | Standard Leaded Glass Mobile Barrier [30 x 24 Leaded Window 33 x 75 Overall Size] |
| LB-3048 | Standard Leaded Glass Mobile Barrier [30 x 48 Leaded Window 33 x 75 Overall Size] |
| LB-3060 | Standard Leaded Glass Mobile Barrier [30 x 60 Leaded Window 33 x 75 Overall Size] |
| LB-3036 | Standard Leaded Glass Mobile Barrier [36 x 30 Leaded Window 33 x 75 Overall Size] |

COMPACT WINDOW MOBILE BARRIERS


LB-1012-2472

LB-1012-3072

LB-1012-3672

| SKU | Part Description |
|--------------|---|
| LB-1012-2472 | Compact Window Leaded Glass Mobile Barrier [12 x 10 Leaded Window 24 x 72 Overall Size] |
| LB-1012-3072 | Compact Window Leaded Glass Mobile Barrier [12 x 10 Leaded Window 30 x 72 Overall Size] |
| LB-1012-3672 | Compact Window Leaded Glass Mobile Barrier [12 x 10 Leaded Window 36 x 72 Overall Size] |

SPECIALTY/ADJUSTABLE MOBILE BARRIERS


LB-3648


LB-3672

LB-7429


| SKU | Part Description | SKU | Part Description |
|---------|--|---------|---|
| LB-3648 | Specialty Leaded Glass Mobile Barrier [36 x 48 Leaded Window 52.5 x 75 Overall Size] | LB-7429 | Adjustable Height Leaded Glass Mobile Barrier [24 x 24 Leaded Window 31 x 41 Opaque Panel 31 x 70 Overall Size] |
| LB-3672 | Specialty Leaded Glass Mobile Barrier [36 x 72 Leaded Window 75.5 x 75 Overall Size] | | |




APRON PEG RACKS




| SKU | Part Description |
|-----------|---|
| ARP-02-61 | Wall Mounted Storage Rack - 2 Pegs 7.5 x 7.5 x 6 [White] |
| ARP-06-61 | Wall Mounted Storage Rack - 6 Pegs 28.5 x 7.5 x 6 [White] |
| ARP-08-61 | Wall Mounted Storage Rack - 8 Pegs 28.5 x 7.5 x 6 [White] |

APRON HANGING RACKS



| SKU | Part Description |
|---------|---|
| AGR-105 | Wall Mounted Apron Hanging Rack - 29 x 7 x 6.25 (L x W x H) [White] |
| AGR-110 | Wall Mounted Apron Hanging Rack - 29 x 7 x 18 (L x W x H) [White] |

SWING-ARM APRON RACKS



| SKU | Part Description |
|--------|---|
| ARW-L5 | Wall Mounted Swing-Arm Apron Rack - 5 Arm (Swing Left) [White] |
| ARW-R5 | Wall Mounted Swing-Arm Apron Rack - 5 Arm (Swing Right) [White] |
| ARW-10 | Wall Mounted Swing-Arm Apron Rack - 5 Arm (Swing L & R) [White] |



MOBILE SWING-ARM APRON RACKS



| SKU | Part Description |
|---------|---|
| ARM-L5 | Mobile Swing-Arm Apron Rack - 5-Arm (Swing Left) [White] |
| ARM-R5 | Mobile Swing-Arm Apron Rack - 5-Arm (Swing Right) [White] |
| ARM-100 | Mobile Swing-Arm Apron Rack - 10-Arm (Swing L & R) [White] |
| ARM-G60 | Mobile Swing-Arm Apron Rack - 10-Arm (Swing L & R) w/ 6 Glove Holders [White] |

Z-STORAGE RACK



| SKU | Part Description |
|----------|--|
| H1978BL | Mobile Z-Base Apron Storage Rack [Black] |
| ARM-UZ10 | Mobile Z-Base Apron Storage Rack w/ 10 Steel Apron Hangers [Black] |

STEEL APRON HANGERS


| SKU | Part Description |
|--------|-------------------------------------|
| ARH-01 | Steel Apron Hanger (Single) |
| ARH-10 | Steel Apron Hanger (Quantity of 10) |



APRON MANAGEMENT


APRON TAGS

| SKU | Part Description |
|--------|---|
| AHDT22 | Apparel 2022 Date Tags - Gray (Pack of 50) |
| AHDT23 | Apparel 2023 Date Tags - Pink (Pack of 50) |
| AHDT24 | Apparel 2024 Date Tags - White (Pack of 50) |
| 500249 | Apron-Hub QR Serial Number Labels (Single Roll of 100 Labels) |
| AQRT | ApronHub Adhesive QR Tag |



FLEXIBLE LEAD BLOCKERS

| SKU | Part Description |
|--------|---|
| FLB-58 | Flexible Lead Blocker - 5 x 8 (1mm Pb Equivalency) |
| FLB-52 | Flexible Lead Blocker - 5 x 12 (1mm Pb Equivalency) |
| FLB-61 | Flexible Lead Blocker - 6 x 10 (1mm Pb Equivalency) |
| FLB-71 | Flexible Lead Blocker - 7 x 11 (1mm Pb Equivalency) |
| FLB-77 | Flexible Lead Blocker - 7 x 17 (1mm Pb Equivalency) |
| FLB-78 | Flexible Lead Blocker Bundle (1 Each of 5 Sizes) |



APRON RECYCLING

Techno-Aide's unique partnership with world renowned transitional company Veolia®, allows us to offer a one-of-a-kind, on-demand, lead apron recycling program. This revolutionary program provides health care facilities of virtually any size, a quick and cost-effective way to dispose of lead and/or lead-equivalent garments properly and safely!

Visit **Techno-Aide.com** for more information!

| SKU | Part Description |
|-------|---|
| RPAK1 | Protective Apparel Recycling Box - (1 Single Box) |
| RPAK5 | Protective Apparel Recycling Box - (5 Boxes) |



Light on Weight. Heavy on Protection.

TechnoAide

Techno-Aide is proud to preset our premium selection of lead and lead equivalent core materials. Our apparel is lined with a Bi-Layer core material that has been proven to provide the highest degree of radiation protection while offering the lightest fitting option for healthcare professionals.

| MATERIAL | REGULAR LEAD | LIGHTWEIGHT LEAD | SUPER-LITE® LEAD-FREE |
|---------------|---|---|---|
| PROTECTION | 0.25mm Pb equivalence: (1-ply) 0.25mm | 0.25mm Pb equivalence: (1-ply) 0.25mm | 0.25mm Pb equivalence: (1-ply) 0.25mm |
| | 0.35mm Pb equivalence: (2-ply) 0.175mm | 0.35mm Pb equivalence: (2-ply) 0.175mm | 0.35mm Pb equivalence: (2-ply) 0.175mm |
| | 0.5mm Pb equivalence: (2-ply) 0.25mm | 0.5mm Pb equivalence: (2-ply) 0.25mm | 0.5mm Pb equivalence: (2-ply) 0.25mm |
| TRANSMISSIONS | -60 kV 0.25mm: 2.5% 0.35mm: 0.9% 0.5mm: 0.32% | -60 kV 0.25mm: 2.44% 0.35mm: 0.82% 0.5mm: 0.23% | -60 kV 0.25mm: 2.39% 0.35mm: 0.86% 0.5mm: 0.19% |
| | -90 kV 0.25mm: 8.75% 0.35mm: 4.95% 0.5mm: 2.48% | -90 kV 0.25mm: 8.72% 0.35mm: 4.93% 0.5mm: 2.49% | -90 kV 0.25mm: 8.74% 0.35mm: 4.91% 0.5mm: 2.45% |
| | -110 kV 0.25mm: 12.69% 0.35mm: 7.67% 0.5mm: 4.39% | -110 kV 0.25mm: 12.67% 0.35mm: 7.65% 0.5mm: 4.38% | -110 kV 0.25mm: 12.60% 0.35mm: 7.63% 0.5mm: 4.36% |
| ATTENUATION | -60 kV 0.25mm: 97.5% 0.35mm: 99.1% 0.5mm: 99.68% | -60 kV 0.25mm: 97.56% 0.35mm: 99.18% 0.5mm: 99.77% | -60 kV 0.25mm: 97.61% 0.35mm: 99.14% 0.5mm: 99.81% |
| | -90 kV 0.25mm: 91.25% 0.35mm: 95.05% 0.5mm: 97.52% | -90 kV 0.25mm: 91.28% 0.35mm: 95.07% 0.5mm: 97.51% | -90 kV 0.25mm: 91.26% 0.35mm: 95.09% 0.5mm: 97.55% |
| | -110 kV 0.25mm: 87.31% 0.35mm: 92.33% 0.5mm: 95.61% | -110 kV 0.25mm: 87.33% 0.35mm: 92.35% 0.5mm: 95.62% | -110 kV 0.25mm: 87.40% 0.35mm: 92.37% 0.5mm: 95.64% |
| TOLERANCES | -10%/+10% | -7%/+2% | -7%/+2% |



Apparel Embroidery

Techno-Aide is proud to offer **FREE** embroidery/monogramming services for all protective apparel products.

Aprons include the following embroidery/monogramming options:

Embroidered Logo on Pockets (*One-Time Digitizing Fee May Apply*)

Up to 2 Lines of Embroidery - 26 Characters Per Line

EMBROIDERY STYLES

STANDARD BLOCK

Sample

ALL CAPS

SAMPLE

SCRIPT

Sample

EMBROIDERY THREAD OPTIONS:




EMBROIDERY LOCATIONS


APRONS/VESTS



SKIRTS



THYROIDS



EMBROIDERY PATCH



ABCDEFGHIJKLMNOPQRSTUVWXYZ
abcdefghijklmnopqrstuvwxyz

Patches include the following embroidery options:

Up to 2 Lines of Embroidery - 26 Characters Per Line

Styles: Standard Block, All Caps



STANDARD SOLID

Lightweight, liquid repellent, filament reinforced fabric



11 - Royal Blue



12 - Navy



20 - Green



21 - Teal



22 - Forest



23 - Gray



54 - Berry



33 - Black



41 - Sunshine



43 - Orange



51 - Red



52 - Hot Pink



53 - Purple



19 - Lime
Green



10 - Sky Blue



55 - Burgundy

NEW COLORS

COMFORT SOLID

Ultra-smooth, lightweight, liquid repellent, filament reinforced fabric



91 - Stratos




93 - Midnight

REINFORCED RIPSTOP


Med. weight, liquid repellent, deluxe filament reinforced fabric




61 - Sapphire



62 - Ruby



63 - Emerald



64 - Onyx

BINDING COLORS

Lightweight, liquid repellent, durable binding



K - Black




G - Gray

PRINTS & PATTERNS


Medium weight, liquid repellent, deluxe filament patterned fabric



69 - Skulls



71 - Leopard



72 - Lot-O-Dots



73 -
Sleepy Sheep



77 - Life Vines



78 - Waterfall



76 - Blue Paws



80 - Red Paws




81 - Black Paws



82 - Pink Camo




85 -
Digital Camo



92 - Ocean



67 - Chevron
Rainbow



68 - Candy
Skulls



70 - Snow
Panther



79 - Shady Plaid




83 - Patriotic
Stars




86 - Night
Camo



87 - Daisies




88 - Splatter



89 - Red Flames



94 - Beach
Party



96 - Funky
Hearts



97 - Trix



98 - Lilac Clouds

NEW DESIGNS

SHOPPING WITH TECHNO-AIDE IS EASY AS 1-2-3

1. Choose the perfect style & design.
2. Enter your favorite fabric color code.
3. Enter your favorite binding color code.



SIZE EXTERIOR COLOR CODE BINDING CODE TETHERED or NON-TETHERED
 WM-5R-1133-K-V-T
 DESIGN MEDIA TYPE INTERIOR COLOR CODE CLOSURE

EXAMPLE - EZ Reverse Vest - Medium | 0.5 Pb EQ | Regular Lead | 11-Royal Blue, 33-Black | Black Binding | Hook & Loop | Thyroid Tether



Size Charts

Techno-Aide's one-of-a-kind **EZ Fit Sizing Tool** quickly and accurately fits people of all shapes and sizes into the perfect Techno-Aide apron size without the need for measuring tapes. However, use the provided measuring charts to verify your fit if a measuring tape is the only available option.


Techno-Aide standard aprons (front, full wrap, vests, & skirts) are fit using two measurements:

- **Width** measured at the **widest point of the body** applicable to apron style
- **Length** measured applicable to apron style

Visit Techno-Aide.com/ez-fit-sizing-tool for step-by-step measuring instructions and to order your sizing tool today!


WIDTHS

FRONT, FULL WRAP, VESTS, & SKIRTS




LENGTHS


APRONS (FRONT & FULL WRAP)



VESTS (STANDARD & REVERSE)



SKIRTS





LAP GUARD HALF APRONS

| Measurement | Petite | Small | Medium | Large | X-Large | 2XL |
|----------------------------|--------|-------|--------|-------|---------|-------|
| Width (at widest point) | 10" | 12" | 20.5" | 24" | 24" | 36" |
| Length (top to bottom) | 10" | 12" | 16" | 17" | 18.5" | 19.5" |

BLANKET GUARDS

| Measurement | Small | Medium | Large | X-Large |
|---------------------------|-------|--------|-------|---------|
| Width (side to side) | 18" | 24" | 24" | 24" |
| Length (top to bottom) | 18" | 24" | 36" | 48" |

PREGNANCY FRONT APRONS

| Measurement | Standard | Full |
|----------------------------|----------|------|
| Width (at chest) | 30" | 34" |
| LENGTH (from shoulders) | 43" | 45" |

PREGNANCY HALF APRONS

| Measurement | Standard |
|------------------------|----------|
| Width (at hips) | 27" |
| Length (from waist) | 27" |

GONADAL GUARDS

| Measurement | Small | Medium | Large |
|---------------------------|-------|--------|-------|
| Width (side to side) | 4" | 6" | 7.5" |
| Length (top to bottom) | 6" | 7" | 9" |
| Strap Length | 8-61" | 8-61" | 8-61" |

BREAST GUARDS

| Measurement | Child | Adult |
|---------------------------|-------|-------|
| Width (side to side) | 5" | 6" |
| Length (top to bottom) | 6" | 7" |
| Strap Length | 8-36" | 8-36" |

THYROID GUARDS

| Measurement | Standard | Mammo | Pediatric |
|----------------------------|----------|-------|-----------|
| Length (strap to strap) | 24" | 24" | 16" |
| Bib Height | 5.75" | 4.75" | 4" |

PEDIATRIC FRONT APRONS

| Measurement | Standard |
|----------------------------|----------|
| Width (at waist) | 18" |
| Length (from shoulders) | 25.5" |

DENTAL DRAPES

| Measurement | Child | Adult |
|--------------------------------|-------|-------|
| Width (at widest point) | 18" | 23.5" |
| Length (from neckline) | 23.5" | 28" |
| Back Flaps (from shoulders) | 7.5" | 11.5" |

PANORAMIC DENTAL DRAPES

| Measurement | Adult |
|----------------------------------|-------|
| Front Width (at widest point) | 24" |
| Front Length (from neckline) | 26" |
| Back Width (at widest point) | 25" |
| Back Length (from neckline) | 15.5" |

SCOLIO GUARDS

| Measurement | Standard |
|---------------------------|----------|
| Width (side to side) | 8" |
| Length (top to bottom) | 26" |
| Neck Strap | 8-36" |
| Waist Strap | 49" |

SIZING TOOL VIDEO

