



A BRIEF GUIDE TO POSITIONING SPONGES

Non-Coated

\$



- Basic (core) piece of porous (open cell) foam, no cover
- No cleanable surface
- Not resistant to fluids
- Not considered durable
- Highly radiolucent
- Artifact free
- Most inexpensive option on average
- Some shapes available in “Stealth” design to reduce the chance of artifacts



Actual Image Under Low Exposure X-Ray

Closed Cell

\$\$



- High density, crosslink (closed cell) foam
- **Cleanable surface; No coating/cover required!**
- Entirely resistant to fluids/liquids
- High durability
- Highly radiolucent
- **Artifact free**
- **Most Cost-Effective option**
- “Stealth” design not necessary




Actual Image Under Low Exposure X-Ray

UNDERSTANDING “STEALTH” SPONGES


NON-STEALTH/ “TRADITIONAL” SPONGES

- Typically feature right angles; most common style for positioning sponges
- Right angles show as artifacts during imaging when a form of coating or covering is applied to core sponge





Coated \$\$\$



- Open cell foam covered with proprietary medical-grade coating
- Cleanable surface
- Resistant to fluids
- More durable than basic (core) foam
- Semi radiolucent
- Can be artifact free
- Second most expensive option on average
- Some shapes available in “Stealth” design to reduce the chance of artifacts



Actual Image Under Low Exposure X-Ray

Vinyl Covered \$\$\$\$




- Open cell foam covered with completely sewn vinyl cover
- Cleanable surface
- Resistant to fluids
- Extreme durability and comfort
- Semi radiolucent
- Not artifact free (sewn seams will show during most imaging)
- Most expensive option on average
- Not available in “Stealth” design



Actual Image Under Low Exposure X-Ray

STEALTH/“FLARED” SPONGES

- Flared sides help significantly reduce chance of artifacts during imaging
- Increase overall size of sponge, typically at base of product





YXBB



YXBL



YXAM



YXAA



YXCA



YXBE



YXAY



YXBU



YXAJ



YXNZ



YXUZ



YXAK



YXWB



CLOSED CELL SPONGES

YXI6	45° Wedge Sponge [30" x 9" x 9"]
YXAA	Finger Block Sponge [5" x 5" x 3"]
YXAD	Adult Bi-Lateral Finger Block Set (L & R) [YXFR & YXFL]
YXAJ	Knee Wedge Sponge [22" x 13.5" x 7"]
YXAK	Circular Head Sponge [Diameter: 7" x 1.5"]
YXAM	Flat Bottom V-Trough/AP & PA Adult Head Holder Sponge [15" x 14" x 4.625"]
YXCM	Large Flat Bottom V-Trough [36" x 14" x 4.625"]
YXAY	Specialty Angular Bolster Sponge [14" x 22.5" x 8"]
YXBB	Medium 45° Spinal Body Wedge Sponge [14" x 7" x 7"]
YXBC	45° Wedge Sponge [7" x 7" x 7"]
YXBD	Large 45° Spinal Body Wedge Sponge [21" x 7" x 7"]
YXBE	34° Wedge Sponge [7" x 4" x 3"]
YXBF	27° Wedge Sponge [10" x 10" x 6"]
YXBG	21° Wedge Sponge [10" x 10" x 4.625"]
YXBI	15° Wedge Sponge [10" x 7" x 2"]
YXBJ	12° Wedge Sponge [9.5" x 7.25" x 2.25"]
YXBL	Rectangle Sponge [10" x 7.5" x 3"]

YXBN	Rectangle Sponge [18" x 14" x 3"]
YXBR	Rectangle Sponge [6" x 4" x 1"]
YXBT	Rectangle Sponge [6" x 4" x 2"]
YXBU	30° x 60° x 90° Multiangle Wedge [29.5" x 10" x 6"]
YXBV	Square Block Sponge [6" x 6" x 2"]
YXFL	Left Lateral Finger Block Sponge [6" x 6.125" x 3.25"]
YXFR	Right Lateral Finger Block Sponge [6" x 6.125" x 3.25"]
YXGF	Rectangle Sponge [12" x 10" x 2"]
YXMB	Rectangle Sponge [10" x 8" x 2"]
YXNZ	Podiatry Axial/Sesamoid Weight Bearing Sponge [13" x 8" x 4"]
YXUZ	Ultrasound Wedge Sponge [20" x 20" x 6"]
YXWB	Weight Bearing Cassette Holder [Overall size: 17.5" x 11.5" x 4" Channel: 1.1875" x 11.5" x 3"]
YXCA	S Extremity Sponge [13" x 7" x 3.5"]
YXCB	M Extremity Sponge [19.5" x 11.5" x 5"]
YXCC	L Extremity Sponge [30.5" x 11" x 7"]
YXCD	XL Extremity Sponge [44" x 14" x 9"]

CLOSED CELL BUNDLES


YB-X-A04	Clinic Sponge Bundle A [YXAA, YXAK, YXBB, YXBI]
YB-X-A05	Clinic Sponge Bundle B [YXAA, YXAK, YXBC, YXBD, YXBE, YXBI, YXBJ, YXBL]
YB-X-A14	General Sponge Bundle 1 [YXAA, YXAK, YXBD, YXBI, YXBN, YXBV, YXMB]
YB-X-A07	General Sponge Bundle A [YXAK, YXBC, YXBE, YXBI, YXBL]
YB-X-A08	General Sponge Bundle B [YXAA, YXAK, YXBC, YXBE, YXBI, YXBJ, YXBL]
YB-X-A13	General Sponge Bundle P [1 ea. - YXAA, YXAK, YXBD, YXBI]
YB-X-A01	Veterinary Sponge Bundle [1 ea. - YXBI, YXBT, YXBR, YXAM, YXBD]
YB-X-A02	Urgent Care Sponge Bundle [1 ea. - YXAA, YXAD, YXAK, YXBC, YXBD, YXBE, YXBJ, YXBL, YXBN]
YB-X-A03	Alternate Urgent Care Sponge Bundle [1 ea. - YXAA, YXAD, YXAK, YXBC, YXBD, YXBE, YXBI, YXBL, YXBN]




Non-Coated



YFBB




YFBO



YFAQ




YFBE




YFAY



YFAK



YFAW



YFAM



YFAB



NON-COATED SPONGES

YBBA	(S) 45° Bariatric Spinal Wedge Sponge [27.5" x 10" x 7"]
YFAA	(S) Adult Oblique Finger Block Sponge [10.5" x 6" x 3"]
YFAB	(S) Bilateral Adult Oblique Finger Block Sponge [12.5" x 10.5" x 3"]
YFAC	(S) Adult Lateral Finger Block Sponge [7.75" x 7.7" x 2"]
YFAD	(S) Adult Lateral Finger Block Sponge [13" x 8" x 3"]
YFAE	(S) Pediatric Lateral Finger Block Sponge [7" x 4" x 1.75"]
YFAK	(S) Circular Head Block [Diameter: 7" (top) x 9" (base) x 1.5"]
YFAL	(S) Head & Neck Support Sponge [16" x 14.5" x 6.5"]
YFAM	(S) AP/PA Adult Head Holder Sponge [15" x 14" x 6.75"]
YFAN	(S) AP/PA Child Head Holder Sponge [11" x 9.75" x 5"]
YFAQ	(S) Child Head Immobilizer Sponge [11.5" x 6.75" x 5.25"]
YFAR	(S) Cervical Head Rest Sponge [16" x 7.5" x 4.5"]
YFAW	Adult Myelogram Block Sponge [22" x 10" x 14"]
YFAX	Child Myelogram Block Sponge [14" x 8" x 10"]
YFAY	Specialty Angular Bolster Sponge [22.5" x 14" x 7.25"]
YFBB	(S) 45° Spinal Wedge Sponge [20.25" x 10" x 7"]
YFBC	(S) 45° Spinal Body Wedge Sponge [13.5" x 10" x 7"]
YFBD	(S) 45° Spinal Body Wedge Sponge [27.75" x 10" x 7"]
YFBE	(S) 34° Wedge Sponge Set [9.75" x 5.5" x 3"]
YFBF	(S) 27° Wedge Sponge [15.5" x 12.75" x 6"]
YFBG	(S) 21° Wedge Sponge [14.25" x 12.25" x 4.5"]


(S) - Stealth Sponges

YFBH	(S) 16° Wedge Sponge Set [11" x 8.5" x 3.5"]
YFBI	(S) 15° Wedge Sponge Set [10" x 12.5" x 3"]
YFBJ	(S) 12° Sponge Set [9.25" x 10.5" x 2.25"]
YFBL	(S) Rectangle Sponge [12.5" x 10.25" x 3"]
YFBN	(S) Rectangle Sponge [20.5" x 16.25" x 3"]
YFBO	(S) Square Sponge [13.5" x 13.5" x 4"]
YFBR	Rectangle Sponge [6" x 4" x 1"]
YFBT	Rectangle Sponge [6" x 4" x 2"]
YFBU	(S) 30° x 60° x 90° Multiangle Wedge Sponge [29.5" x 10" x 6"]
YFBW	CT Headrest Sponge [12.5" x 11" x 10"]
YFBY	(S) Rectangle Sponge [27.5" x 21.5" x 4"]
YFCA	(S) Small Extremity Sponge [12.5" x 7" x 3"]
YFCB	(S) Medium Extremity Sponge [19.75" x 11.5" x 5"]
YFCC	(S) Large Extremity Sponge [24" x 11" x 7"]
YFCD	(S) X-Large Extremity Sponge [36" x 14" x 9"]
YFDB	30° Cardiac Wedge Sponge [27" x 28" x 16"]
YFEA	(S) 45° Spinal Wedge Sponge [24" x 13.75" x 9.5"]
YFEE	(S) 45° Echocardiography Wedge Sponge [38.25" x 13" x 9"]
YFEM	Rectangle Sponge [24" x 15" x 2.5"]
YFFD	45° Spinal Body Wedge Sponge [21" x 7" x 7"]


NON-COATED BUNDLES

YB-F-A01	Thin Rectangular Pad (Set of 2) [24" x 18" x 0.25"]
YB-F-A02	Urgent Care Sponge Bundle [1 ea. - YFAA, YFAD, YFAK, YFBC, YFBD, YFBE, YFBJ, YFBL, YFBN]
YB-F-A04	Clinic Sponge Bundle A [1 ea. - YFBB, YFAA, YFAK, YFBI]
YB-F-A05	Clinic Sponge Bundle B [1 ea. - YFAA, YFAK, YFBC, YFBE, YFBI, YFBJ, YFBI, YFBD]
YB-F-A06	CT Positioning Sponge Bundle [1 ea. - YFBB, YFBC, YFBE, YFBE, YFBE, YFBI, YFBJ, YFBI, YFBD]
YB-F-A07	General Sponge Bundle A [1 ea. - YFAK, YFBC, YFBE, YFBI, YFBL]
YB-F-A08	General Sponge Bundle B [1 ea. - YFAA, YFAK, YFBC, YFBE, YFBI, YFBJ, YFBL]

YB-F-A10	General Sponge Bundle D [1 ea. - YFAL, YFAN, YFAR, YFBE, YFBJ, YFBL]
YB-F-A12	General Sponge Bundle G [1 ea. - YFAA, YFAK, YFBD, YFBE, YFBI, YFBL]
YB-F-A13	General Sponge Bundle P [1 ea. - YFAA, YFAK, YFBD, YFBI]
YB-F-A16	Ortho Sponge Bundle A [1 ea. - YFAA, YFAM, YFAR, YFBD, YFBL]
YB-F-A19	Pediatric Sponge Bundle [1 ea. - YFAC, YFAE, YFAK, YFAN, YFAQ, YFAX, YFBT]
YB-F-A20	Wedge Sponge Bundle A [1 ea. - YFBC, YFBF, YFBG, YFBI]
YB-F-A21	Wedge Sponge Bundle B [1 ea. - YFBD, YFBE, YFBE, YFBG]




YCBA




YCBL




YCBF



YCEP




YCCA



YCAV



YCAT




YCAG



YCAA



YCMG



YCAW



YCAQ



YCAK



YMCP



COATED SPONGES

YCAA	(S) Adult Oblique Finger Block Sponge [10.5" x 6" x 3"]
YCAC	(S) Bilateral Adult Oblique Finger Block Sponge [12.5" x 10.5" x 3"]
YCAC	(S) Bilateral Pediatric Oblique Finger Block Sponge [7.75" x 7.7" x 2"]
YCAD	(S) Adult Lateral Finger Block Sponge [13" x 8" x 3"]
YCAE	(S) Pediatric Lateral Finger Block Sponge [7" x 4" x 1.75"]
YCAG	(S) Bilateral Leg Immobilizer Sponge [21" x 20" x 6"]
YCAH	(S) Knee Arthrogram Block Sponge [19" x 14.5" x 8"]
YCAJ	(S) 26" Knee Wedge Support Sponge [21" x 19" x 8"]
YCAK	(S) Circular Head Block [Diameter: 7" (top) x 9" (base) x 1.5"]
YCAL	(S) Head & Neck Support Sponge [14" x 16" x 6.5"]
YCAM	(S) AP/PA Adult Head Holder Sponge [15" x 14" x 6.75"]
YCAN	(S) AP/PA Child Head Holder Sponge [11" x 9.75" x 5"]
YCAP	(S) Adult Head Immobilizer Sponge [15.25" x 9" x 7"]
YCAQ	(S) Child Head Immobilizer Sponge [11.5" x 6.75" x 5.25"]
YCAR	(S) Cervical Head Rest Sponge [16" x 7.5" x 4.5"]
YCAT	(S) Bariatric Torso Sponge [16.5" x 18.5" x 7"]
YCAW	Adult Myelogram Block Sponge [22" x 10" x 14"]
YCAX	Child Myelogram Block Sponge [14" x 8" x 10"]
YCAZ	(S) Specialty Angular Bolster Sponge [22.5" x 14" x 7.25"]
YCBA	(S) 45° Bariatric Spinal Wedge Sponge [27.5" x 10" x 7"]
YCBB	(S) 45° Spinal Wedge Sponge [20.25" x 10" x 7"]
YCBC	(S) 45° Spinal Body Wedge Sponge [13.5" x 10" x 7"]
YCBD	(S) 45° Spinal Body Wedge Sponge [27.75" x 10" x 7"]
YCBE	(S) 34° Wedge Sponge Set [9.75" x 5.5" x 3"]
YCBF	(S) 27° Wedge Sponge [15.5" x 12.75" x 6"]
YCBG	(S) 21° Wedge Sponge [14.25" x 12.25" x 4.5"]
YCBH	(S) 16° Wedge Sponge Set [11" x 8.5" x 3.5"]
YCBI	(S) 15° Wedge Sponge Set [10" x 12.5" x 3"]
YCBJ	(S) 12° Wedge Sponge Set [9.25" x 10.5" x 2.25"]
YCBL	(S) Rectangle Sponge [12.5" x 10.25" x 3"]
YCBM	(S) Bariatric Rectangle Sponge [14.5" x 10.5" x 3"]
YCBN	(S) Rectangle Sponge [20.5" x 16.25" x 3"]
YCBO	(S) Square Sponge [13.5" x 13.5" x 4"]
YCBP	(S) Bariatric Decubitus Rectangle Sponge [27.5" x 21.5" x 4"]
YCBQ	Rectangle Sponge [6" x 4" x 0.5"]

(S) - Stealth Sponges | *non-returnable. More info upon request.

YCBR	Rectangle Sponge [6" x 4" x 1"]
YCBS	Rectangle Sponge [6" x 4" x 1.5"]
YCBT	Rectangle Sponge [6" x 4" x 2"]
YCBU	(S) 30° x 60° x 90° Multiangle Wedge [29.5" x 10" x 6"]
YCBV	Square Sponge [6" x 6" x 2"]
YCBW	CT Headrest Sponge [12.5" x 11" x 10"]*
YCBY	(S) Rectangle Sponge [27.5" x 21.5" x 4"]
YCCA	(S) Small Extremity Sponge [12.5" x 7" x 3"]
YCCB	(S) Medium Extremity Sponge [19.75" x 11.5" x 5"]
YCCC	(S) Large Extremity Sponge [30.5" x 11" x 7"]
YCCD	(S) X-Large Extremity Sponge [44.5" x 14" x 9"]
YCCG	45° Bariatric Body Wedge Sponge [30" x 9" x 9"]
YCDA	45° Cardiac Wedge Sponge [28" x 19" x 14"]
YCDB	30° Cardiac Wedge Sponge [27" x 28" x 16"]
YCEA	(S) 45° Spinal Wedge Sponge [33" x 13.75" x 9.5"]
YCEB	(S) 45° Wedge Sponge [30" x 10" x 7.5"]
YCEE	(S) 45° XL Echocardiography Wedge Sponge [38.25" x 13" x 9"]
YCEG	65° Wedge Sponge [25" x 16" x 8.25"]
YCEM	Rectangle Sponge [24" x 15" x 2.5"]
YCEP	Bariatric Endo-Ultrasound Wedge Sponge [16" x 16" x 4.5"]
YCEQ	Perineal Endovaginal/Transvaginal Ultrasound Wedge Sponge [20" x 20" x 3"]
YCET	Occipital Head Support Sponge [13" x 12.75" x 4.5"]*
YCEU	30° x 60° x 90° Multiangle Bariatric Wedge Sponge [29.5" x 10" x 6"]
YCEV	CT/PET Head Holder Sponge Set [13.5" x 14" x 6.75"]*
YCEW	Philips CT Headrest Sponge [12" x 12.5" x 6.5"]*
YCFD	45° Spinal Wedge Sponge [21" x 7" x 7"]
YCFL	(S) Left Lateral Finger Sponge [5-6.5" x 7.5" x 3t"]
YCFR	(S) Right Lateral Finger Sponge [5-6.5" x 7.5" x 3"]
YCFQ	Decubitus Block & Image Receptor Plate/Cassette Holder [Lower: 36" x 16" x 4 Back 36" x 10" x 11" Channel: 36" x 0.5" x 3"]
YCGG	Rectangle Sponge [18" x 8" x 1.5"]
YCMG	Contoured Torso Sponge [18.5" x 17" x 5"]
YCUZ	Endo-Ultrasound Wedge Sponge [20" x 20" x 6"]
YMCS	(S) Small Leg Rest w/ Knee Cutouts [18" x 14.5" x 4"]
YMCM	(S) Medium Leg Rest w/ Knee Cutouts [16.25" x 12.5" x 7"]
YMCP	(S) Large Leg Rest w/ Knee Cutouts [24" x 18" x 10"]



YVFO



YVFJ



YVFQ



YVFA



YVFB



YVWS



YVAM



YB-V-A01



COATED BUNDLES

YCAS	Knee Arthroplasty Stress Assembly Sponge Set
YB-C-A22	Bariatric Sponge Bundle [1ea. - YCAT, YCBA, YCBM, YCBP]
YB-C-A04	Clinic Sponge Bundle A [1 ea. - YCAA, YCAK, YCBB, YCBI]
YB-C-A05	Clinic Sponge Bundle B [1 ea. - YCAA, YCAK, YCBC, YCBD, YCBE, YCBI, YCBJ, YCBL]
YB-C-A06	CT Sponge Bundle A [1 ea. - YCBB, YCBC, YCBE, YCBF, YCBG, YCBI, YCBJ, YCBW]
YB-C-A07	General Sponge Bundle A [1 ea. - YCAK, YCBC, YCBE, YCBI, YCBL]
YB-C-A08	General Sponge Bundle B [1 ea. - YCAA, YCAK, YCBC, YCBE, YCBI, YCBJ, YCBL]
YB-C-A09	General Sponge Bundle C [1 ea. - YCAA, YCAM, YCAR, YCAY, YCBL, YCAH]
YB-C-A10	General Sponge Bundle D [1 ea. - YCAL, YCAN, YCAR, YCBF, YCBJ, YCBL]
YB-C-A11	General Sponge Bundle E [1 ea. - YCAA, YCAK, YCAM, YCBB, YCBG, YCBI, YCBN, YCBU]
YB-C-A12	General Sponge Bundle G [1 ea. - YCAA, YCAK, YCBD, YCBE, YCBI, YCBL]
YB-C-A13	General Sponge Bundle P [1 ea. - YCAA, YCAK, YCBD, YCBI]

YB-C-A16	Orthopedic Sponge Bundle A [1ea. - YCAA, YCAM, YCAR, YCBD, YCBL]
YB-C-A17	Orthopedic Sponge Bundle B [1 ea. - YCAL, YCAN, YCAR, YCBC, YCBG, YCBJ, YCBL]
YB-C-A18	Orthopedic Sponge Bundle D [1 ea. - YCAA, YCAR, YCBD, YCBI, YCBL]
YB-C-A19	Pediatric Sponge Bundle [1 ea. - YCAC, YCAE, YCAK, YCAN, YCAQ, YCAX]
YB-C-A23	Rectangle Sponge Bundle [3 ea. - YCBQ, YCBR, YCBS, YCBT]
YB-C-A24	Scanning Support Sponge Bundle [10 ea. - YCBT]
YB-C-A25	General Sponge Bundle D [1 ea. - YCAY, YCEE, YCUZ]
YB-C-A26	Ultrasound Breast or Paracentesis Sponge Bundle [1 ea. - YCBG, YCFD, YCBU]
YB-C-A15	General Sponge Bundle 1 [2 ea. - YCAK ; 1 ea. - YCBD, YCBI, YCBL, YCBT]
YB-C-A20	Wedge Sponge Bundle A [1 ea. - YCBC, YCBF, YCBG, YCBI]
YB-C-A21	Wedge Sponge Bundle B [1 ea. - YCBD, YCBE, YCBF, YCBG]

VINYL SPONGES

YVAM	Flat Bottom V-Trough/AP & PA Adult Head Holder Sponge [15" x 14" x 6.75"]
YVBI	15" Wedge Sponge Set [10" x 7" x 2"]
YVBR	Rectangle Sponge [6" x 4" x 1"]
YVBT	Rectangle Sponge [6" x 4" x 2"]
YVDA	45° Cardiac Wedge Sponge [28" x 14" x 14"]
YVF40	Knee Wedge Bolster Sponge [17" x 13" x 6"]
YVFA	Cylinder Bolster Sponge [24" x 8" x 8"]
YVFB	Semi-Cylinder Bolster Sponge [18" x 8" x 6"]
YVFC	20° Wedge Bolster Sponge [24" x 18" x 6"]
YVFD	45° Wedge Bolster Sponge [16" x 8" x 8"]
YVFE	20° Incline Wedge Bolster Sponge [32" x 20" x 12"]
YVFF	50° Incline Wedge Bolster Sponge [32" x 20" x 12"]
YVFG	Rectangle Bolster Sponge [12" x 8" x 6"]

YVFH	Rectangle Bolster Sponge [14" x 12" x 3"]
YVFJ	Rectangle Bolster Sponge [12" x 8" x 3"]
YVFK	Rectangle Bolster Sponge [20" x 13" x 3"]
YVFM	Knee Wedge Bolster Sponge [25" x 16" x 8.25"]
YVFO	45° Spinal Wedge Bolster Sponge [21" x 7" x 7"]
YVFP	Rectangle Bolster Sponge [24" x 18" x 4"]
YVFQ	DEXA Leg Bolster Sponge [20" x 20" x 12"]
YVFU	30° x 60° x 90° Wedge Sponge [24" x 10" x 6"]
YVFX	Comfort Pillow Body Wedge [25" x 16" x 8.25"]
YVFZ	Endo-Ultrasound Wedge Sponge [20" x 20" x 6"]
YVWA	Specialty Stretcher & Wheelchair Wedge Sponge [24.5" x 14" x 5.5-14"]
YVWS	Specialty Wheelchair Bolster Sponge [17" x 13" x 6"]

VINYL BUNDLE

YB-V-A01	Vinyl Covered Veterinary Sponge Bundle (1 Each - YVBI, YVBT, YVBR, YVAM, YVFO)
----------	---



Sandbags, Sandsocks, & MRI Positoners



**HEAVY-GAUGE VINYL SQUARE
SANDBAGS:
STANDARD HANDLES**



**HEAVY-GAUGE VINYL SQUARE
SANDBAGS:
AC JOINT HANDLES**



FLEXIBLE VINYL SANDSOCKS



SANDSOCK BUNDLES



**HEAVY-GAUGE VINYL RECTANGULAR
SANDBAGS: STANDARD HANDLES**



**HEAVY-GAUGE VINYL RECTANGULAR
SANDBAGS: AC JOINT HANDLES**



**MRI RECTANGLE POSITIONER SET
GP0132**



**3lb MRI SANDBAG
SPG1003-01**



**MRI SANDBAG POSITIONER BUNDLE
SPG1013-01**

NYLON MRI POSITIONERS

GP0132	MRI Rectangle Positioner Set (2) [14" x 3" x 1"]
GP0152	MRI Rectangle Positioner Set (2) [5" x 3" x 1"]
GP0162	MRI Wedge Positioner Set (2) [10" x 10" x 3"]
GP0172	MRI Wedge Positioner Set (2) [14" x 4" x 1.5"]
GP0182	MRI Rectangle Positioner Set (2) [4" x 3" x 2"]
GP0192	MRI Curved Wedge Positioner Set (2) [14" x 4" x 1.5"]
GP0202	MRI Semi-Round Neck Roll Positioner [7" x 3" x 2.3"]
GE1211	MRI Knee Wedge Positioner [28.2" x 16" x 8"]
GP0102	MRI Positioner Bundle [1 ea. - MFP-01 - 06]

MRI VINYL SANDBAGS

SBG1003-01	Square Sandbag: 3lb (8 x 8) Standard Handle [Black]
SBG1005-01	Rectangle Sandbag: 5lb (5 x 15) Standard Handle [Black]
SBG1007-01	Rectangle Sandbag: 7lb (9 x 12) Standard Handle [Black]
SBG1010-01	Rectangle Sandbag: 10lb (9 x 14) Standard Handle [Black]
SBG1012-01	Square Sandbag: 12lb (12 x 12) Standard Handle [Black]
SBG1013-01	Sandbag Bundle: One of each size [Black]



SANDSOCKS

1lb (2 x 8)	Standard	8lb (7 x 22)	Standard
Black	SS-01R-BLK-NH	Black	SS-08R-BLK-NH
Tan	SS-01R-TAN-NH	Tan	SS-08R-TAN-NH
3lb (4 x 12)	Standard	10lb (6 x 20)	Standard
Black	SS-03R-BLK-NH	Black	SS-10R-BLK-NH
Tan	SS-03R-TAN-NH	Tan	SS-10R-TAN-NH
4lb (8 x 10)	Standard	11lb (6 x 28)	Standard
Black	SS-04R-BLK-NH	Black	SS-11R-BLK-NH
Tan	SS-04R-TAN-NH	Tan	SS-11R-TAN-NH
5lb (6 x 15)	Standard	12lb (6 x 25)	Standard
Black	SS-05R-BLK-NH	Black	SS-12R-BLK-NH
Tan	SS-05R-TAN-NH	Tan	SS-12R-TAN-NH
7lb (6 x 18)	Standard	15lb (6 x 40)	Standard
Black	SS-07R-BLK-NH	Black	SS-15R-BLK-NH
Tan	SS-07R-TAN-NH	Tan	SS-15R-TAN-NH

SANDSOCK BUNDLES

BUNDLE 1	Standard	BUNDLE 3	Standard
Black	SSK-01-BLK-NH	Black	SSK-03-BLK-NH
Tan	SSK-01-TAN-NH	Tan	SSK-03-TAN-NH
BUNDLE 2	Standard	BUNDLE 4	Standard
Black	SSK-02-BLK-NH	Black	SSK-04-BLK-NH
Tan	SSK-02-TAN-NH	Tan	SSK-04-TAN-NH
BUNDLE 2	Standard		
Black	SSK-02-BLK-NH		
Tan	SSK-02-TAN-NH		

RECTANGLE SANDBAGS

1lb (4 x 5)	Standard	ACJoint	7lb (6 x 14)	Standard	ACJoint
Black	SB-01R-BLK-ST	SB-01R-BLK-AC	Black	SB-07R-BLK-ST	SB-07R-BLK-AC
Blue	SB-01R-BLU-ST	SB-01R-BLU-AC	Blue	SB-07R-BLU-ST	SB-07R-BLU-AC
2lb (4 x 9)	Standard	ACJoint	10lb (7 x 16)	Standard	ACJoint
Black	SB-02R-BLK-ST	SB-02R-BLK-AC	Black	SB-10R-BLK-ST	SB-10R-BLK-AC
Blue	SB-02R-BLU-ST	SB-02R-BLU-AC	Blue	SB-10R-BLU-ST	SB-10R-BLU-AC
3lb (5 x 9)	Standard	ACJoint	12lb (8 x 18)	Standard	ACJoint
Black	SB-03R-BLK-ST	SB-03R-BLK-AC	Black	SB-12R-BLK-ST	SB-12R-BLK-AC
Blue	SB-03R-BLU-ST	SB-03R-BLU-AC	Blue	SB-12R-BLU-ST	SB-12R-BLU-AC
5lb (6 x 12)	Standard	ACJoint			
Black	SB-05R-BLK-ST	SB-05R-BLK-AC			
Blue	SB-05R-BLU-ST	SB-05R-BLU-AC			

SQUARE SANDBAGS

5lb (9 x 9)	Standard	ACJoint
Black	SB-05S-BLK-ST	SB-05S-BLK-AC
Blue	SB-05S-BLU-ST	SB-05S-BLU-AC
7lb (10 x 10)	Standard	ACJoint
Black	SB-07S-BLK-ST	SB-07S-BLK-AC
Blue	SB-07S-BLU-ST	SB-07S-BLU-AC
10lb (11 x 11)	Standard	ACJoint
Black	SB-10S-BLK-ST	SB-10S-BLK-AC
Blue	SB-10S-BLU-ST	SB-10S-BLU-AC

SANDBAG BUNDLES

Bundle 1	Standard	Pediatric	Standard
Black	SBK-01-BLK-ST	Black	SBK-03-BLK-ST
Blue	SBK-01-BLU-ST	Blue	SBK-03-BLU-ST
Bundle 2	Standard	Femoral/Angio	Standard
Black	SBK-02-BLK-ST	Black	SBK-04-BLK-ST
Blue	SBK-02-BLU-ST	Blue	SBK-04-BLU-ST

SANDBAG PAIRS

5lb Pair (9 x 9)	Standard	ACJoint
Black	SBPR-05S-BLK-ST	SBPR-05S-BLK-AC
Blue	SBPR-05S-BLU-ST	SBPR-05S-BLU-AC
10lb Pair (11 x 11)	Standard	ACJoint
Black	SBPR-10S-BLK-ST	SBPR-10S-BLK-AC
Blue	SBPR-10S-BLU-ST	SBPR-10S-BLU-AC



**TWO STEP WEIGHT BEARING PLATFORM
(COMFORT STEPS)
WMP-16**



**THREE STEP WEIGHT BEARING PLATFORM
(REINFORCED STEPS)
WMP-28**



**PATIENT FOOT IMMOBILIZER
FIL-12**



**ADJUSTABLE LEG SUPPORT
LGS**



**DUAL PLATFORM BARIATRIC STEP STOOL
31220**

SPECIALTY POSITIONING AIDS

FIL-12	Patient Foot Immobilizer
LGS	Adjustable Leg Support - Bariatric Base [15" x 18"]
LGS-S	Adjustable Leg Support - Standard Base [9" x 12"]
RHI-AB	Patient Head Immobilizer [Coated Sponge - 6 x 6 x 2]
RIS-60-BLK	Heavy-Gauge Vinyl Patient Immobilizer Strap [60" x 10"]
RIS-72-BLK	Heavy-Gauge Vinyl Patient Immobilizer Strap [72" x 10"]
RIS-84-BLK	Heavy-Gauge Vinyl Patient Immobilizer Strap [84" x 10"]
RIS-96-BLK	Heavy-Gauge Vinyl Patient Immobilizer Strap [96" x 10"]
RLC-36	Lateral Chest Image Support Positioner [36"]
RLC-48	Lateral Chest Image Support Positioner [48"]

STANDARD & MRI PLATFORM TOOLS

T-50	Standard Single Platform Step Stool with Handrail
T-6150	Single Platform Bariatric Step Stool with Handrail
31220	Dual Platform Bariatric Step Stool with Handrail
FSA3001-01	MRI Single Platform Step Stool with Handrail
FSM726-01	MRI Dual Platform Step Stool with Handrail
WMP-16	Two Step Weight Bearing Imaging Platform (Closed Cell)
WMP-22	Two Step Weight Bearing Imaging Platform (Plastic)
WMP-14	Three Step Weight Bearing Imaging Platform (Closed Cell)
WMP-28	Three Step Weight Bearing Imaging Platform (Plastic)