

# Improve Workflow

## Upgrade To Digital with AcquiDR



### Cost effective retro-fit DR

- Digital Flat Panel Detector (FPD)
- Imaging Workstation with Software
- Portable/Wireless for Maximum Flexibility
- Easy Interface with Analogue X-ray Systems

Improve workflow. Take high quality images.  
Use a simple user interface – all when you  
partner with Imperial Imaging to upgrade  
to AcquiDR series flat panel detectors.



**Upgrade  
Today**

Be confident that your imaging  
services are right for your business  
today—and in the future.

**678.733.9990 | [Imperialimaging.com](http://Imperialimaging.com)**

# Improve Performance

## Choose AcquiDR from Imperial Imaging

The Imperial Advantage is more than just a product purchase. Our service makes you confident that your equipment empowers you to provide better diagnosis and patient care.

Imperial shares expertise in all facets of imaging products. Compare brands. Meet compliance. Receive training and worry-free installation. Secure financing. Imperial manages it all!



### Fully Integrated DR System

- Full integration with X-ray generator
- Full featured imaging software
- Excellent digital image quality
- DICOM 3.0 Compatible

## Great for Retrofit and New DR Systems

- Wireless – with power cable included for fixed panel options
- Stable iSync and Automatic Exposure Detection (AED)
- Dual band (2.4 and 5 GHz) wireless support with easy sharing
- Long battery life and smart workflow
- Direct deposition CSI for better image quality with lower dose

### Radmax Software Features:

- Automatic Masking
- Study List Filter
- History Search
- Image Zoom
- Grid Suppression
- DICOM Store
- PMS Integration
- Image Merge
- Veterinary specific tools (VHS, TPLO)
- Anatomical Programming (APR)
- Edge Enhancement Processing

### Specifications

Model Name	MANO 4336W	MANO 4343W	MANO 4343T
Design			
Detector Type	Wireless	Wireless	Wired
Image Size	17x14 inches (43x36cm)	17x17 inches (43x43cm)	17x17 inches (43x43cm)
Scintillator	CsI		
Detector Technology	Amorphous Silicon		
AD Conversion	16 Bit		
Digital Software	Radmax		



**Upgrade  
Today**

Imperial helps you select equipment that meets your needs today and for the future.

**678.733.9990**  
**Imperialimaging.com**